Metropolitan Medical Response System (MMRS)

Equipment Technical Data and Information

1999 (revised)

MMRS Functions

- Medical treatment/Management
- Hospital/Public Health
- Incident Command Liaison
- Patient Transportation
- Agent Detection
- Decontamination
- Pharmacology
- Law Enforcement Intelligence/Security
- Mental Health Services
- Non-Survivor Disposition
- Training

MMRS Equipment

- Medical Emergency
- Personal Protective Equipment
- Agent Detection
- Patient Decontamination
- Communications
- Pharmacology Supplies
- Training Materials
- Transportation*

^{*}Correlated to Community Requirement

The Research Planning, Inc. (RPI) Metropolitan Medical Response System (MMRS) Development Support Team continues to expand and strengthen its association with equipment manufacturers and vendors, offering equipment or services that may be useful in an MMRS response. It is our pleasure to pass this information on to the MMRS Development cities as part of our effort to continue "Sharing the Knowledge."

The MMRS Development Support Team has assembled the "Equipment Technical Data and Information" book to assist in the MMRS Development and Equipment Identification process. The information in this book has been provided by manufacturers or vendors, and the use of trade or manufacturers' names does not constitute a recommendation or endorsement of any commercial product. Any equipment specifications listed are subject to change without notice. The MMRS Support Team respectfully reminds all readers of this book that it is the user's responsibility to determine the suitability or fitness of the items for a particular use or application.

NOTICE

Since the conditions of use of any of the equipment identified in this book are beyond the control of the U.S. Department of Health and Human Services (DHHS) and RPI, DHHS and RPI make no warranties, expressed or implied, regarding the suitability or fitness of any of the identified equipment for a particular use and assume no obligation or liability in connection with any use of this information.

Total Washington Metropolitan Area (WMA) Metropolitan Medical Strike Team (MMST) Equipment Requirements

| | WMA MMST | | | | |
|----------------------------------------------------------------------------------------|----------------|-------------|-----------------------------------------|--|--|
| | Unit Price | Quantit | Qty Total | | |
| Developed Brode of the English and (DDE) | | у | | | |
| Personal Protective Equipment (PPE) | | | | | |
| Section 1 – Level A ENTRY Team (PPE) | | | | | |
| KAPPLER RESPONDER PLUS Level A Suits | \$ 543.0 | 8 | \$ 4,344.0 | | |
| | 0 | | 0 | | |
| Bio-Pak 240 Rebreather (4-hr SCBA) | 4,695.00 | 8 | 37,560.00 | | |
| Pairs of North Butyl Gloves Pairs North Silver Shields | 13.45 4.90 | 10 10 | 134.50 49.00 | | |
| Pairs Best "N" Dex Nitrile Gloves | 1.80 | 10 | 18.00 | | |
| ILC Dover Cool Vest | 189.00 | 10 | 1,890.00 | | |
| 120 20101 0001 1001 | 100.00 | Total | \$ 43,995.50 | | |
| | | | | | |
| Section 2 – Level C Decon Team (PPE) | _ | | | | |
| Powered Air Purifying Respirator (PAPR) | \$ 425.7 | 129 | \$ 54,921.75 | | |
| Niced Pottorice (Eytro) | 5 109.20 | 40 | 4 COE CO | | |
| Nicad Batteries (Extra) 10 Bank PAPR Battery Charger | 596.54 | 43 7 | 4,695.60 4,175.78 | | |
| TYVE-F Decon Suits (Color Gray) | 81.00 | 129 | 10,449.00 | | |
| Pairs of North Butyl Gloves | 13.45 | 129 | 1,735.05 | | |
| Pairs North Silver Shields | 4.90 | 129 | 632.10 | | |
| 129 pr Best "N" Dex Nitrile Gloves (25 Pr/Box) | 14.50 | 7 | 101.50 | | |
| Pairs of HAZMAX 16" Kneeboots | 43.30 | 129 | 5,585.70 | | |
| TOPPS Nomex IIIA Coveralls | 82.94 | 129 | 10,699.26 | | |
| Racal Operations and Recognition Vests | 99.99 | 129 | 12,898.71 | | |
| Advent Safety Helmet | 60.00 | 129 | 7,740.00 | | |
| Safety Glasses Hazardous Materials (HAZMAT) Gear Bag | 3.35 39.99 | 129 129 | 432.15 5,158.71 | | |
| Booklets of M-8 Detection Paper (2 Per Person) | 150.00 | 5 | 750.00 | | |
| Dooklets of W o Detection 1 aper (2 1 cm of cison) | 150.00 | Total | \$119,975.31 | | |
| | | | *************************************** | | |
| Section 3 – Replacement PPE Stored in Cache | | | | | |
| TYVEK "F" Decon Suit (3 Extra Suits Per Team Position) | \$ 75.0 | 125 | \$ 9,375.0 | | |
| TVA/EIX "E" Dance Out (O.E. v. O. v. D. T. D. vv.) | 0 | | 0 | | |
| TYVEK "F" Decon Suit (3 Extra Suits Per Team Position) TYVEK "F" Decon Suit (Training) | 81.00 75.00 | 4 43 | 324.00 | | |
| Pairs of North Butyl Gloves | 13.45 | 129 | 3,225.00 1,735.05 | | |
| Pairs North Silver Shields | 4.90 | 258 | 1,264.20 | | |
| 258 pr Best "N" Dex Nitrile Gloves (25 Pr/Box) | 14.50 | 12 | 174.00 | | |
| Spare Oxygen Bottles for Bio Pak 240 | 500.00 | 5 | 2,500.00 | | |
| Kit to Service Bio Paks | 600.00 | 1 | 600.00 | | |
| Lime Pak Refill in 44 gal Drum (11 Refills of Bio Pak) | 252.50 | 2 | 505.00 | | |
| Replacement Ice Pack for Cool Vest | 75.00 | 1 | 75.00 | | |
| Boxes GMAH Org Vapor/HEPA Carts (6 Carts/Box–3 Carts/PAPR) | 90.84 | 22 Total | 1,998.48 | | |
| | | Total | \$ 21,775.73 | | |
| PPE Total (Sections 1, 2, and 3 combined) | | | \$185,746.54 | | |

Total WMA MMST Equipment Requirements (Continued)

| | | WM. | A MMST |
|----------------------------------------------------------------|---------------|--------------|-----------------|
| | Unit Price | Quantit y | Qty Total |
| Communications Equipment | | , | |
| MTS2000 Model II Portable 800 MHz Analog Battery, Antenna | \$ 1,972.00 | 40 | \$ 78,880.0 |
| | | | 0 |
| MTS2000 Model III Portable Analog Battery, Antenna | 2,227.00 | 3 | 6,681.00 |
| Swivel Leather Carrying Case | 34.00 | 43 | 1,462.00 |
| Extra Ultra-High Capacity Batteries | 88.00 | 43 | 3,784.00 |
| Throat Mikes for Radios | 515.00 | 8 | 4,120.00 |
| Multi-Unit (Six Pockets) Rapid Charger | 573.75 | 6 | 3,442.50 |
| Single Unit Rapid Charger | 106.25 | 1 | 106.25 |
| Spectra W7 Mobile 800 MHz 35W /Digital | 2,807.00 | 1 | 2,807.00 |
| Spectra C7 Mobile 403-433 MHz 40W Analog | 1,717.00 | 1 | 1,717.00 |
| Spectra C7 Mobile 450-482 MHz 40W Analog | 1,717.00 | 1 | 1,717.00 |
| Spectra C7 Mobile 482-512 MHz 40W Analog | 1,717.00 | 1 | 1,717.00 |
| Spectra C7 Mobile 136-174 MHz 50W Analog | 1,630.00 | 1 | 1,630.00 |
| Maxtrac 300 Mobile 29.7-36 MHz 60W Analog | 750.00 | 1 | 750.00 |
| Maxtrac 300 Mobile 42-50 MHz 60W Analog | 750.00 | 1 | 750.00 |
| Program Portables & Mobiles, Provide & Install Mobile/NIDA Box | 7,737.00 | 1 | 7,737.00 |
| Pager Service for Team Members (Annual) | 119.00 | 0 | |
| Tagor convict roam monitors (rumaal) | 110.00 | Total | \$117,300.75 |
| | | | |
| <u>Detection Equipment</u> | | | |
| Box M-8 Detection Paper (50 Booklets/Box) | \$ 150.00 | 1 | \$ 150.00 |
| Rolls M-9 Detection Paper | 35.00 | 2 | 70.00 |
| M-256 Detection Kit | 120.00 | 1 | 120.00 |
| M-18 Detection Kit | 310.00 | 2 | 620.00 |
| Draeger Kit | 852.75 | 1 | 852.75 |
| Extra Draeger Tubes | 367.55 | 1 | 367.55 |
| Deluxe Pump Kit for Draeger Kit | 374.50 | 1 | 374.50 |
| Surface Acoustic Wave (SAW) MiniCAD | 5,495.00 | 1 | 5,495.00 |
| Chemical Agent Monitor (CAM) w/Alarms and Tips | 7,500.00 | 2 | 15,000.00 |
| CAM Audible Alarm | 150.00 | 2 | 300.00 |
| Ludlum Digital Scaler Ratemeter 2241 | 675.00 | 1 | 675.00 |
| Ludlum Digital Scaler Ratemeter 2241-2 | 725.00 | 1 | 725.00 |
| 44-6 Side Window GM Probe (Beta/Gamma) | 130.00 | 1 | 130.00 |
| 44-6 Side Window GM Probe (Beta/Gamma) | 145.00 | 1 | 145.00 |
| 44-9 Pancake GM Detector (Alpha) | 175.00 | 1 | 175.00 |
| Self-Reading Dosimeters | 110.00 | 10 | 1,100.00 |
| Dosimeter Charger Thormolyminoscopt Dosimeter (TLD) Type 1/1/9 | 95.00 6.30 | 1 43 | 95.00 270.90 |
| Thermoluminescent Dosimeter (TLD) Type W8 | 6.30 | _ | |
| | | Total | \$ 26,665.70 |

Total WMA MMST Equipment Requirements (Continued)

| | | WM | A MMST |
|----------------------------------------------------|------------|---------|--------------|
| | Unit Price | Quantit | Qty Total |
| | 0111111100 | у | aty rotal |
| Decontamination Equipment | | , | |
| M-295 Individual Decon Kit | \$ 28.00 | 60 | \$ 1,680.00 |
| M-295 Individual Decon Kit | 41.50 | 69 | 2,863.50 |
| Decontamination System – Primary +\$2,000 Shipping | 65,000.00 | 1 | 65,000.00 |
| Decontamination System Supplies – Secondary | 2,053.06 | 1 | 2,053.06 |
| Decon Shower | 520.00 | 1 | 520.00 |
| Di-chlorite Powder (28 lb Buckets) | 31.95 | 2 | 63.90 |
| Chlorine Measurement Kit | 52.30 | 1 | 52.30 |
| 95 Gal Poly Over Pak Drums | 200.95 | 5 | 1,004.75 |
| Cartons Belongings Bags | 149.80 | 3 | 449.40 |
| Paper Hospital Gown | 1.20 | 400 | 480.00 |
| Cartons Cadaver Bags | 295.20 | 3 | 885.60 |
| Cartons Cadaver Dags | 255.20 | Total | \$ 75,052.51 |
| | | Total | Ψ 73,032.31 |
| Medical Equipment | | | |
| SharpSafe Containers (100 Per Case) | \$ 310.00 | 2 | \$ 620.00 |
| Adult 18 GA. | 1.14 | 1000 | 1,140.00 |
| Adult 20 GA. | 1.14 | 1000 | 1,140.00 |
| PED. 22 GA. | 1.14 | 500 | 570.00 |
| Bacteriostatic Fluid | 3.25 | 100 | 325.00 |
| Butterfly 21 GA. | 0.50 | 500 | 250.00 |
| Lactated Ringers | 2.80 | 1000 | 2,800.00 |
| 18 GA. Needle | 7.95 | 10 | 79.50 |
| 20 GA. Needle | 7.95 | 10 | 79.50 |
| Sphygmomanometer (C) | 22.85 | 15 | 342.75 |
| Sphygmomanometer (I) | 22.85 | 10 | 228.50 |
| Sphygmomanometer (A) | 14.00 | 20 | 280.00 |
| E.T. (2.5) w/o Cuff | 2.25 | 20 | 45.00 |
| E.T. (3.5) w/o Cuff | 2.25 | 20 | 45.00 |
| NASO Airways 6mm | 4.70 | 20 | 94.00 |
| NASO Airways 8mm | 4.70 | 20 | 94.00 |
| Surg Mask w/Face Shield | 27.50 | 8 | 220.00 |
| O.B. Kit | 8.25 | 5 | 41.25 |
| Syringes, 1cc | 0.14 | 1000 | 140.00 |
| Syringes, 3cc | 0.15 | 1000 | 150.00 |
| Syringes, 5cc | 0.20 | 1000 | 200.00 |
| Syringes, 20cc | 0.65 | 1000 | 650.00 |
| Syringes, 30cc | 0.80 | 1000 | 800.00 |
| Syringes, 60cc | 1.25 | 1000 | 1,250.00 |
| Boxes Needles (21 g., 1"), 100 per Box | 7.95 | 30 | 238.50 |
| Oxygen Manifold – Minilator | 199.00 | 8 | 1,592.00 |
| Ambu Bags | 19.95 | 20 | 399.00 |
| O2 Tubing 7' | 0.46 | 50 | 23.00 |
| Nasal Cannula 7' | 1.50 | 50 | 75.00 |
| Non-Rebreather Mask | 1.40 | 50 | 70.00 |
| Carry Case | 96.50 | 3 | 289.50 |
| Stethoscope | 7.25 | 20 | 145.00 |
| ' | | Total | \$14,416.50 |

Total WMA MMST Equipment Requirements (Continued)

| | | WM | A MMST |
|-------------------------------------------------|-------------|------------|--------------|
| | Unit Price | Quantit | Qty Total |
| | | У | |
| Ancillary Equipment | | | |
| Computer, Fax Modem, and Cabin and Med Software | \$ 4,225.00 | 2 | \$ 8,450.00 |
| Miscellaneous General Support Items | 2,700.00 | 1 | 2,700.00 |
| Portable Tent/Shelter | 6,519.25 | 1 | 6,519.25 |
| 20-Foot Utility Trailer (Equip/Command) | 4,608.01 | 1 | 4,608.01 |
| Diesel Generator | 3,359.00 | 1 | 3,359.00 |
| Deployment Packaging | 1,939.94 | 1 | 1,939.94 |
| Personnel Command Vests | 13.00 | 15 | 195.00 |
| Boards, Phones, etc. | 1,000.00 | 1 | 1,000.00 |
| Color Coded Tarps | 20.83 | 6 | 124.98 |
| Raven Stretchers | 346.75 | 10 | 3,467.50 |
| Rolls Barrier Tape | 13.90 | 10 | 139.00 |
| | | Total | \$ 32,502.68 |
| | | | |
| | ship | ping/taxes | \$ 7,419.73 |
| | TOTAL | | \$459,104.41 |

Prices are estimates and do not include handling or administrative costs. Figures are to be used as a planning tool only.

MMST Equipment for 50 Persons

| | | Total | |
|------------------------------------------|--------------|--------------|--------------|
| | Unit Price | Requirements | Qty Total |
| MMST INDIVIDUAL EQUIPMENT | | • | |
| | | | |
| Level C PPE | | | |
| Individual Personal Protection Equipment | | | |
| RRPAS with NiCAD | \$ 514.39.00 | 50 | \$ 25,719.50 |
| NiCAD Batteries (Extra) | 109.20 | 50 | 5,460.00 |
| 10 Bank PAPR Battery Charger | 596.54 | 5 | 2,982.70 |
| TYVEK-F Decon Suits (Color Gray) | 81.00 | 100 | 8,100.00 |
| Pairs of North Butyl Gloves | 13.45 | 100 | 1,345.00 |
| Pairs North Silver Shields | 4.90 | 100 | 490.00 |
| Best "N" Dex Nitrile Gloves (25 Pr/Box) | 14.50 | 4 | 58.00 |
| Pairs of HAZMAX 16" Kneeboots | 43.30 | 50 | 2,165.00 |
| TOPPS Nomex IIIA Coveralls | 82.94 | 50 | 4,147.00 |
| "KIWI" Safety Helmet | 125.00 | 50 | 6,250.00 |
| Safety Glasses | 3.35 | 50 | 167.50 |
| HAZMAT Gear Bag | 39.99 | 50 | 1,999.50 |
| Boxes GMAH Org Vapor/HEPA Carts | 90.84 | 25 | 2,271.00 |
| (6 Carts/Box – 3 Carts/PAPR) | | | |
| Individual PPE Equipment | | | \$ 61,155.20 |
| Individual PPE Equipment Per Person | | | \$ 1,223.10 |

| | Uı | nit Price | Total Requirements | (| Qty Total |
|---------------------------------------------------------------|----|-----------|-----------------------|----|-----------|
| MMST TEAM EQUIPMENT | | | | | |
| Level A ENTRY Team (PPE) KAPPLER RESPONDER PLUS Level A Suits | \$ | 543.00 | 12 | \$ | 6,516.00 |
| Bio-Pak 240 Rebreather (4-hr SCBA) | | 4,695.00 | 12 | | 56,340.00 |
| Pairs of North Butyl Gloves | | 13.45 | 12 | | 161.40 |
| Pairs North Silver Shields | | 4.90 | 12 | | 58.80 |
| Best "N" Dex Nitrile Gloves (25 Pr/Box) | | 1.80 | 1 | | 1.80 |
| ILC Dover Cool Vest | | 189.00 | 12 | | 2,268.00 |
| HAZMAT Gear Bag | | 39.99 | 12 | | 479.88 |
| Spare Oxygen Bottles for Bio Pak 240 | | 500.00 | 12 | | 6,000.00 |
| Kit to Service Bio Paks | | 600.00 | 1 | | 600.00 |
| Lime Pak Refill in 44 gal Drum (11 Refills of Bio Pak) | | 252.50 | 2 | | 505.00 |
| Replacement Ice Pack for Cool Vest | | 75.00 | 1 | | 75.00 |
| Level A Team Equipment | | | | \$ | 73,005.88 |

| | | Total | |
|-------------------------------------------|------------|--------------|--------------|
| | Unit Price | Requirements | Qty Total |
| | | | |
| Detection Equipment | | | |
| Box M-8 Detection Paper (50 Booklets/Box) | \$ 150.00 | 1 | \$ 150.00 |
| Rolls M-9 Detection Paper | 35.00 | 2 | 70.00 |
| M-256 Detection Kit | 120.00 | 1 | 120.00 |
| M-18 Detection Kit | 310.00 | 1 | 310.00 |
| Draeger Kit | 852.75 | 1 | 852.75 |
| Extra Draeger Tubes | 367.55 | 1 | 367.55 |
| Deluxe Pump Kit for Draeger Kit | 374.50 | 1 | 374.50 |
| SAW MiniCAD | 5,495.00 | 1 | 5,495.00 |
| CAM w/Alarms and Tips | 8,000.00 | 1 | 8,000.00 |
| 4-Gas Portable Monitor | 1,700.00 | 1 | 1,700.00 |
| Photoionization Detector | 1,325.00 | 1 | 1,325.00 |
| Ludlum Digital Scaler Ratemeter 2241-2 | 725.00 | 1 | 725.00 |
| 44-6 Side Window GM Probe (Beta/Gamma) | 145.00 | 1 | 145.00 |
| 44-9 Pancake GM Detector (Alpha) | 175.00 | 1 | 175.00 |
| Self-Reading Dosimeters | 110.00 | 12 | 1,320.00 |
| Dosimeter Charger | 95.00 | 1 | 95.00 |
| TLD Type W8 | 6.30 | 50 | 315.00 |
| Detection Equipment Package | | | \$ 21,539.80 |

| | Unit Price | Total Requirements | Qty Total |
|---------------------------------------------|------------|-----------------------|--------------|
| Decontamination Equipment | | | |
| M-295 Individual Decon Kit | \$ 41.50 | 50 | \$ 2,075.00 |
| Decontamination System Supplies – Secondary | 2,053.06 | 1 | 2,053.06 |
| Decon Shower | 520.00 | 1 | 520.00 |
| Sodium Hypochlorite Powder (15-lb Buckets) | 58.12 | 8 | 464.96 |
| 95-gal Poly Over Pak Drums | 200.95 | 5 | 1,004.75 |
| Cartons Belongings Bags | 149.80 | 3 | 449.40 |
| Paper Hospital Gowns | 1.20 | 400 | 480.00 |
| Cartons Cadaver Bags | 295.20 | 4 | 1,180.80 |
| Decontamination Equipment Sub-Package | | | \$ 8,227.97 |
| Decontamination System – Primary | 65,000.00 | 1 | 65,000.00 |
| Transportation and Shipping (Average) | 2,000.00 | 1 | 2,000.00 |
| Decontamination Equipment Package | | | \$ 75,227.97 |

| | Unit F | Price | Total Requirements | Qty Total |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------|--------------------------------------|-----------------------|-----------------------------------------------------------------|
| Ancillary Equipment Miscellaneous General Support Items 20-Foot Utility Trailer (Equip/Command) Raven Stretchers Field Support Packaging Ancillary Equipment Package | 5,0 | 000.00 091.85 065.00 000.00 | 1 1 10 1 | \$ 3,000.00 5,091.85 3,650.00 2,500.00 \$ 14,241.85 |
| Total MMST Team Equipment Total MMST Individual Equipment | | | Shipping | 4,000.00 \$188,015.50 \$ 61,155.20 \$249,170.70 |
| Prices are estimates and do not include handling of | r adminis | | quipment Total costs. | \$249,170.70 |

| PRODUCT | PRODUCT CODE | TEAM MEMBER PERSONAL CACHE | IMMEDIATE ON-SCENE TEAM CACHE | SECONDARY CALL-UP CACHE | TOTAL INVENTORY | UNIT ORDER SIZE (Package Size) | TOTAL UNITS | ESTIMATED UNIT COST | ESTIMATED PRODUCT COST | ANTICIPATED PRODUCT USE (On-Scene or Hospital Resupply) |
|-----------------------------------------------------|------------------|----------------------------|-------------------------------------|----------------------------|--------------------|-----------------------------------|-------------|------------------------|------------------------|---------------------------------------------------------|
| ANALGESICS | | | | | | | | | | |
| etoralac (Toradol) Injection 30mg/mL; 1mL | 6505-01-357-5375 | 0 | 25 | 25 | 50 | 10 | 5 | \$50.99 | \$254.95 | ОН |
| CONTROLLED ANALGESICS | | | | | | | | | | |
| utorphanol Injection 2mg/mL | 6505-00-000-2978 | 0 | 25 | 25 | 50 | 10 | 5 | \$31.28 | \$156.40 | ОН |
| orphine Sulphate Injection 10mg/mL | 6505-00-149-0113 | 0 | 25 | 25 | 50 | 10 | 5 | \$5.57 | \$27.85 | ОН |
| ANESTHETICS | | | | | | | | | | |
| etamine HCl Injection 10mg/mL; 20mL Vial | NA | 0 | 20 | 30 | 50 | 10 | 5 | \$66.46 | \$332.30 | OH |
| idocaine 1% Injection; 20mL Vial | 6505-00-000-2819 | 0 | 20 | 30 | 50 | 25 | 2 | \$12.43 | \$24.86 | ОН |
| idocaine 4% Viscous; 50mL | 6505-00-772-0157 | 0 | 20 | 30 | 50 | 1 | 50 | \$8.53 | \$426.50 | ОН |
| ANTICONVULSANTS | | | | | | | | | | |
| henytoin Injection 50mg/mL; 5mL | 6505-00-139-4348 | 0 | 20 | 30 | 50 | 10 | 5 | \$9.77 | \$48.85 | OH |
| henobarbital Injection 130mg/mL | 6505-00-812-2556 | 0 | 0 | 50 | 50 | 10 | 5 | \$21.87 | \$109.35 | OH |
| TOPICAL ANTI-INFECTIVES | | | | | | | | | | |
| olysporin Ointment; 3gm Unit of Use Packets | NA | 0 | 144 | 144 | 288 | 144 | 2 | \$13.83 | \$27.66 | ОН |
| ilver Sulfadiazine (Silvadene) Cream 1%; 85gm | NA | 0 | 20 | 80 | 100 | 1 | 100 | \$5.73 | \$573.00 | ОН |
| CARDIOVASCULAR PRODUCTS | | | | | | | | | | |
| tenolol 5mcg/ml Injection | NA | 0 | 20 | 30 | 50 | 1 | 50 | \$10.36 | \$518.00 | ОН |
| alcium Chloride Injection 10% 10mL | 6505-00-009-4968 | 0 | 20 | 30 | 50 | 10 | 5 | \$14.01 | \$70.05 | ОН |
| igoxin Injection 0.25mg/mL; 2mL amps | 6505-00-531-7761 | 0 | 20 | 30 | 50 | 10 | 5 | \$13.43 | \$67.15 | ОН |
| opamine 250mg/600ml Injection | NA | 0 | 20 | 30 | 50 | 10 | 5 | \$19.08 | \$95.40 | ОН |
| pinephrine 1:10,000 (0.1mg/mL); 10mL | 6505-01-149-8559 | 0 | 30 | 20 | 50 | 25 | 2 | \$36.58 | \$73.16 | OH |
| eparin 25,000/250ml D5W | NA | 0 | 5 | 7 | 12 | 12 | 1 | \$33.27 | \$33.27 | OH |
| idocaine 1% Injection; 2mL Prefilled Syringes (PFS) | 6505-01-105-4179 | 0 | 20 | 30 | 50 | 10 | 5 | \$45.93 | \$229.65 | ОН |
| lagnesium Sulfate Injection; 1gm | NA | 0 | 20 | 30 | 50 | 10 | 5 | \$63.36 | \$316.80 | ОН |
| itroglycerin Injection 5mg/mL; 10mL | NA | 0 | 5 | 10 | 15 | 10 | 1.5 | \$4.59 | \$6.89 | ОН |
| itroglycerin Tablets 1/150gr (bottles) | 6505-00-687-3663 | 0 | 100 | 100 | 200 | 100 | 2 | \$4.01 | \$8.02 | ОН |
| otassium Chloride Injection; 10meg; 5mL | NA | 0 | 10 | 20 | 30 | 10 | 3 | \$48.12 | \$144.36 | ОН |

| PRODUCT | PRODUCT CODE | TEAM MEMBER PERSONAL CACHE | IMMEDIATE ON- SCENE TEAM CACHE | SECONDARY CALL UP CACHE | TOTAL INVENTORY | UNIT ORDER SIZE | TOTAL UNITS | ESTIMATED UNIT COST | ESTIMATED PRODUCT COST | ANTICIPATED PRODUCT USE (On Scene or Hospital Resupply) |
|--------------------------------------------------|------------------|----------------------------|--------------------------------------|-------------------------|--------------------|-----------------|-------------|------------------------|---------------------------|---------------------------------------------------------|
| odium Bicarbonate Injection 50mEq; 50mL | 6505-00-000-2921 | 0 | 10 | 20 | 30 | 10 | 3 | \$16.28 | \$48.84 | ОН |
| treptokinase 1.5mu Injection | NA | 0 | 0 | 10 | 10 | 1 | 10 | \$184.53 | \$1,845.30 | ОН |
| ELECTROLYTES | | | | | | | | | | |
| atorade (Oral Rehydration Electrolyte Powder) | NA | 0 | 84 | 60 | 144 | 12 | 12 | \$18.94 | \$227.28 | 0 |
| PSYCHOTHERAPEUTIC PRODUCT | | | | | | 10 | 0 | | | |
| aloperidol Injection 5mg/mL | 6505-00-268-8530 | 0 | 20 | 20 | 40 | 10 | 4 | \$3.17 | \$12.68 | ОН |
| PULMONARY PRODUCTS | | | | | | | | | | |
| trovent Inhalers | 6505-01-256-1948 | 0 | 10 | 10 | 20 | 1 | 20 | \$18.09 | \$361.80 | ОН |
| eclomethasone Inhaler | 6505-01-240-0588 | 0 | 10 | 10 | 20 | 1 | 20 | \$14.09 | \$281.80 | ОН |
| piPen (or AnaKits) | 6505-01-043-6795 | 0 | 20 | 20 | 40 | 1 | 40 | \$8.14 | \$325.60 | ОН |
| SERUMS, TOXOIDS, AND VACCINES | | | | | | | | | | |
| amma Globulin; 5mL | NA | 0 | 0 | 10 | 10 | 0 | 0 | | \$0.00 | Н |
| CHEMICAL AND NERVE GAS ANTIDOTES | | | | | | | | | | |
| .tropine/2-PAM Antidote Kit | NA | 378 | 600 | 2022 | 3000 | 1 | 3000 | \$13.01 | \$39,030.00 | ОН |
| ana Kit (Diazepam) | NA | 126 | 200 | 674 | 1000 | 1 | 1000 | \$11.25 | \$11,250.00 | ОН |
| illy Cyanide Antidote Kit | 6505-00-000-2966 | 0 | 20 | 20 | 40 | 1 | 40 | \$106.70 | \$4,268.00 | ОН |
| 291 Skin Decontaminating Kit | NA | 0 | 0 | 1000 | 1000 | 5000 | 0.2 | \$202.80 | \$40.56 | Н |
| ethylprednisolone 80mg/mL; 1mL | 6505-00-000-2890 | 0 | 20 | 20 | 40 | 10 | 4 | \$21.22 | \$84.88 | ОН |
| ralidoxime (2-PAM) | | | 12 | 12 | 24 | 1 | 24 | \$85.00 | \$2,040.00 | ОН |
| BIOLOGICAL ANTIDOTES | | | | | | | | | | |
| mantadine - 200 mg Capsules | 6505-00-148-4624 | 0 | 0 | 400 | 400 | 100 | 4 | \$10.30 | \$41.20 | Н |
| mphotericin B Injection | 6505-01-084-9453 | 0 | 0 | 10 | 10 | 1 | 10 | \$6.78 | \$67.80 | Н |
| icillin Cartridge Needle Units 1ml 600,000 Units | 6505-00-133-4448 | 0 | 0 | 50 | 50 | 10 | 5 | \$54.87 | \$274.35 | Н |
| eftazidine 2gm Vial | NA | 0 | 0 | 25 | 25 | 10 | 2.5 | \$174.58 | \$436.45 | Н |
| iprofloxacin 500mg Capsules | 6505-01-333-4155 | 0 | 0 | 1000 | 1000 | 100 | 10 | \$170.66 | \$1,706.60 | Н |
| oxycycline 200mg Capsules | 6505-01-153-4335 | 0 | 0 | 500 | 500 | 500 | 1 | \$29.00 | \$29.00 | Н |

| PRODUCT | PRODUCT CODE | TEAM MEMBER PERSONAL CACHE | IMMEDIATE ON- SCENE TEAM CACHE | SECONDARY CALL UP CACHE | TOTAL INVENTORY | UNIT ORDER SIZE (Package Size) | TOTAL UNITS | ESTIMATED UNIT COST | ESTIMATED PRODUCT COST | ANTICIPATED PRODUCT USE (On-Scene or Hospital Resupply) |
|-----------------------------------------------------------------------|------------------|----------------------------------|--------------------------------------|-------------------------|--------------------|-----------------------------------|-------------|------------------------|------------------------|---------------------------------------------------------|
| Pentavalent Antitoxin and Toxoid (IND) | NA | 0 | 0 | | 0 | | | \$0.00 | \$0.00 | Н |
| Rifampin 300mg Capsules | 6505-00-165-6575 | 0 | 0 | 500 | 500 | 100 | 5 | \$44.31 | \$221.55 | Н |
| Streptomycin 400mg/ml Injection | NA | 0 | 0 | | 0 | 10 | 0 | \$0.00 | \$0.00 | Н |
| Superactivated Charcoal | NA | 0 | 0 | 5 | 5 | 1 | 5 | \$7.15 | \$35.75 | Н |
| Trimethoprim (TMP) and Sulfmethoxicol (SMX) 5ml Injection | NA | 0 | 0 | 50 | 50 | 10 | 5 | \$4.17 | \$20.85 | Н |
| Vibramycin 100mg Capsules | 6505-00-369-7296 | 0 | 0 | 100 | 100 | 50 | 2 | \$2.97 | \$5.94 | Н |
| Virazole - 6gms Powder / 100ml Vial (available only under IND status) | NA | 0 | 0 | 4 | 4 | 4 | 1 | \$4,031.15 | \$4,031.15 | Н |
| RADIATION SPECIFIC | | | | | | | | | | |
| Potassium Iodide, Tablets | 6505-01-116-8198 | 0 | 0 | 4494 | 4494 | 14 | 321 | \$1.06 | \$340.26 | Н |
| ANCILLARY SUPPLIES | | | | | | | | | | |
| Ambu Bag, Adult | NA | 0 | 50 | 4 | 54 | 6 | 9 | \$176.40 | \$1,587.60 | ОН |
| Ambu Bag, Pediatric | NA | 0 | 25 | 5 | 30 | 6 | 5 | \$255.22 | \$1,276.10 | ОН |
| Bandage Scissors 7" | 6515-00-363-8840 | 0 | 20 | 20 | 40 | 1 | 40 | \$2.19 | \$87.60 | ОН |
| Biohard Bag | NA | 0 | 150 | 350 | 500 | 250 | 2 | \$20.57 | \$41.14 | ОН |
| Bite Blocks | 6515-01-100-8656 | 0 | 50 | 0 | 50 | 10 | 5 | \$4.56 | \$22.80 | ОН |
| Flashlight w/Batteries | 6532-00-781-3671 | 0 | 10 | 20 | 30 | 1 | 30 | \$5.81 | \$174.30 | ОН |
| Gauze 4" x 4" | 6510-00-782-2698 | 0 | 200 | 400 | 600 | 200 | 3 | \$3.48 | \$10.44 | ОН |
| Handi-Wipes | NA | 0 | 2000 | 0 | 2000 | 100 | 20 | \$2.17 | \$43.40 | 0 |
| Heparin Locks 10units/3ml | NA | 0 | 200 | 0 | 200 | 200 | 1 | \$453.60 | \$453.60 | 0 |
| Intravenous Catheter & Needle Unit, 14G x 2.25" | NA | 0 | 100 | 100 | 200 | 200 | 1 | \$359.38 | \$359.38 | OH |
| Intravenous Catheter & Needle Unit, 18G x 2" | NA | 0 | 100 | 50 | 150 | 50 | 3 | \$79.43 | \$238.29 | OH |
| Intravenous Catheter & Needle Unit, 22G x 2" | NA | 0 | 50 | 50 | 100 | 50 | 2 | \$79.54 | \$159.08 | OH |
| Intravenous Set, Butterfly, 12", 21G, 3/4" | NA | 0 | 50 | 70 | 120 | 120 | 1 | \$80.40 | \$80.40 | OH |
| IV Administration Set, 78", w/Clamp, Vented | NA | 0 | 24 | 72 | 96 | 48 | 2 | \$33.12 | \$66.24 | ОН |

| PRODUCT | PRODUCT CODE | TEAM MEMBER PERSONAL CACHE | IMMEDIATE ON- SCENE TEAM CACHE | SECONDARY CALL UP CACHE | TOTAL INVENTORY | UNIT ORDER SIZE (Package Size) | TOTAL UNITS | ESTIMATED UNIT COST | ESTIMATED PRODUCT COST | ANTICIPATED PRODUCT USE (On-Scene or Hospital Resupply) |
|-------------------------------------------------------------------------|------------------|----------------------------|--------------------------------------|-------------------------|--------------------|-----------------------------------|-------------|------------------------|---------------------------|---------------------------------------------------------|
| IV Administration Set, MICRO, w/clamp vented | NA | 0 | 24 | 72 | 96 | 48 | 2 | \$76.32 | \$152.64 | ОН |
| Kling Gauze, 3" | 6510-00-717-5094 | 0 | 10 | 10 | 20 | 1 | 20 | \$0.45 | \$9.00 | OH |
| Laryngoscope (disposable) large (123mm x 25.7mm blade) #3 MacIntosh | 56903 | 0 | 25 | 20 | 45 | 1 | 45 | \$36.48 | \$1,641.60 | ОН |
| Laryngoscope (disposable) small (77.5mm x 18.2mm blade) #2 MacIntosh | 56901 | 0 | 10 | 20 | 30 | 1 | 30 | \$36.48 | \$1,094.40 | ОН |
| Laryngoscope Light Source (handle for C batteries) | 60500 | 0 | 5 | 25 | 30 | 1 | 30 | \$29.75 | \$892.50 | OH |
| Latex Gloves Sterile (Lg) | 6515-01-149-8840 | | 200 | 200 | 400 | 50 | 8 | \$21.96 | \$175.68 | OH |
| Latex Gloves Sterile (Med) | 6515-01-151-1790 | 0 | 200 | 200 | 400 | 50 | 8 | \$15.80 | \$126.40 | ОН |
| Latex Gloves Sterile (Xlg) | 6515-01-149-8843 | 0 | 200 | 200 | 400 | 50 | 8 | \$22.67 | \$181.36 | ОН |
| Mask, Surgical, Disposable, with Face Shield (and eye shield) | NA | 0 | 50 | 50 | 100 | 100 | 1 | \$105.12 | \$105.12 | ОН |
| Nas-Gastric Tube #16 | NA | 0 | 20 | 0 | 20 | 1 | 20 | \$238.80 | \$4,776.00 | |
| Nasal Cannulas | 6515-01-166-8099 | 0 | 20 | 50 | 70 | 1 | 70 | \$0.45 | \$31.50 | OH |
| Nasopharmyngeal Airways, 8.0mm ID, 10.5mm OD | NA | 0 | 5 | 5 | 10 | 10 | 1 | \$57.36 | \$57.36 | OH |
| Nasopharyngeal Airways, 6.0mm ID, 8.0mm OD | NA | 0 | 5 | 5 | 10 | 10 | 1 | \$57.36 | \$57.36 | ОН |
| Needle, Disposable, 18G x 1.5" | 6515-00-754-2834 | 0 | 100 | 100 | 200 | 100 | 2 | \$3.27 | \$6.54 | ОН |
| Needle, Disposable, 22G x 1" | 6515-00-754-2835 | 0 | 100 | 100 | 200 | 100 | 2 | \$2.74 | \$5.48 | OH |
| Needle, Disposable, 22G x 1.5" | 6515-00-059-5782 | 0 | 100 | 100 | 200 | 100 | 2 | \$2.62 | \$5.24 | OH |
| Needle, Disposable, 25G x 5/8" | 6515-00-655-5751 | 0 | 100 | 100 | 200 | 100 | 2 | \$3.64 | \$7.28 | ОН |
| Obstetrical Kit, Emergency | NA | 0 | 1 | 1 | 2 | 1 | 2 | \$72.00 | \$144.00 | ОН |
| Pressure Infuser (500 mL) | AD-500 | 0 | 12 | 0 | 12 | 1 | 12 | \$298.80 | \$3,585.60 | 0 |
| SAM Splints | NA | 0 | 12 | 12 | 24 | 12 | 2 | \$202.32 | \$404.64 | ОН |
| Sodium Chloride 0.9%, 500 ml, Plastic Bags | NA | 0 | 96 | 48 | 144 | 24 | 6 | \$16.42 | \$98.52 | ОН |
| Sphygmomanometer, Set, Child | NA | 0 | 6 | 6 | 12 | 1 | 12 | \$10.50 | \$126.00 | ОН |

| PRODUCT | PRODUCT CODE | TEAM MEMBER PERSONAL CACHE | IMMEDIATE ON- SCENE TEAM CACHE | SECONDARY CALL UP CACHE | TOTAL INVENTORY | UNIT ORDER SIZE (Package Size) | TOTAL UNITS | ESTIMATED UNIT COST | ESTIMATED PRODUCT COST | ANTICIPATED PRODUCT USE (On-Scene or Hospital Resupply) |
|------------------------------------------------|------------------|----------------------------|--------------------------------------|-------------------------|--------------------|-----------------------------------|-------------|------------------------|------------------------|---------------------------------------------------------|
| Sphygmomanometer, Set, Infant | NA | 0 | 6 | 6 | 12 | 1 | 12 | \$11.82 | \$141.84 | ОН |
| Sylets 10FR | NA | 0 | 20 | 0 | 20 | 20 | 1 | \$89.75 | \$89.75 | OH |
| Sylets 16FR | NA | 0 | 20 | 0 | 20 | 20 | 1 | \$89.75 | \$89.75 | OH |
| Syringe/Needle, Disposable, 20 g xl.6", 5-6 ml | 6515-00-985-7183 | 0 | 100 | 100 | 200 | 100 | 2 | \$20.62 | \$41.24 | ОН |
| Syringe/Needle, Disposable, 3 ml, 21G x 1 1/2" | 6515-00-462-7248 | 0 | 100 | 100 | 200 | 100 | 2 | \$4.36 | \$8.72 | ОН |
| Tape 1" | 6510-01-060-6369 | 0 | 12 | 24 | 36 | 12 | 3 | \$16.63 | \$49.89 | ОН |
| Tape 2" | 6510-01-095-9285 | 0 | 12 | 24 | 36 | 6 | 6 | \$16.57 | \$99.42 | ОН |
| Tape 3" | 6510-01-060-6370 | 0 | 12 | 24 | 36 | 4 | 9 | \$16.60 | \$149.40 | OH |
| Triage Tags | NA | 0 | 500 | 0 | 500 | 50 | 10 | \$12.81 | \$128.10 | OH |
| Tube, Endotracheal with cuff, 5.0mm | NA | 0 | 5 | 5 | 10 | 10 | 1 | \$38.82 | \$38.82 | OH |
| Tube, Endotracheal with cuff, 6.5mm | NA | 0 | 5 | 5 | 10 | 10 | 1 | \$38.82 | \$38.82 | OH |
| Tube, Endotracheal with cuff, 7.5mm | NA | 0 | 5 | 5 | 10 | 10 | 1 | \$38.82 | \$38.82 | OH |
| Tube, Endotracheal with cuff, 8.0mm | NA | 0 | 5 | 5 | 10 | 10 | 1 | \$38.82 | \$38.82 | ОН |
| Tube, Endotracheal without Cuff, 2.0mm | NA | 0 | 5 | 5 | 10 | 10 | 1 | \$27.00 | \$27.00 | ОН |
| Tube, Endotracheal without cuff, 3.0mm | NA | 0 | 5 | 5 | 10 | 10 | 1 | \$27.00 | \$27.00 | ОН |
| Tube, Endotracheal without cuff, 4.0mm | NA | 0 | 5 | 5 | 10 | 10 | 1 | \$27.00 | \$27.00 | ОН |
| TOTAL COST FOR THE MMST PHARMACEUTICALS | | | | | | | | | \$89,891.09 | |

PERSONAL PROTECTIVE EQUIPMENT

PERSONAL PROTECTIVE EQUIPMENT: Respirators

Equipment Type:

Negative Pressure, Full Facepiece Respirators

Make / Model:

3M "6000" Series and "7000" Series

Description:

These are tight-fitting, negative pressure respirators that utilize dual bayonet-mount filter cartridges. The lightweight, silicone rubber facepiece features a one-piece polycarbonate lens and a separate nose cup to reduce internal lens fogging.

Features / Benefits

- Large single lens for wide field of vision
- Easy to don and doff
- 6000 Series Adaptable to dual airline system, fully elastic 4-point head harness
- 7000 Series Adaptable to Powered Air Purifying Respirator (PAPR) or supplied air systems, fully elastic six-point head harness, integral speaking diaphragm

Specifications

Three sizes: Small, Medium, and Large

Weight

- 6000 Series 14.5 oz.
- 7000 Series 20.5 oz.

Accessories (available at extra cost)

- Spectacle kit
- Fit Test Kit with Bayonet-type Attachment Adapter
- Face Shield Cover
- Tinted Lens Cover
- 6000-series mask available as a component of the "Chem / Bio Response-Pak," available through TriCon Environmental, Inc. (see below)

Approvals

- National Institute of Occupational Safety and Health (NIOSH) approvals:
 - Filter part # 2091 P-100 particulate efficiency
 - Filter part # 60921 Organic Vapor cartridge with P-100 particulate filter
 - Filter part # 60923 Organic vapor, chlorine, hydrogen chloride and sulfur dioxide (escape only) or hydrogen fluoride with P-100 particulate filter
 - Filter part # 60926 Organic vapor, chlorine, hydrogen chloride, chlorine dioxide, sulfur dioxide, hydrogen sulfide (escape only), ammonia/methylamine, or hydrogen fluoride cartridge with P-100 particulate filter
 - Lens meets requirements of American National Standards Institute (ANSI) Z87.1 1989 (Face and Eye Protection)

Manufacturer / Distributor

3M Occupational Health and Safety Division St. Paul, MN 55133 TriCon Environmental, Inc. Huntsville, AL

(800) 854-4334 / www.tricon-env.com

(800) 896-4223 / www.3m.com



3M "6000" Series and "7000" Series Negative Pressure, Full Facepiece Respirators

Equipment Type:

Negative Pressure, Full Facepiece Respirators

Make / Model:

3M "SGE-400"

Description:

This is a tight-fitting, negative pressure respirator utilizing a single standard North Atlantic Treaty Organization (NATO)-thread pattern filter cartridge. The lightweight, silicone rubber facepiece features a one-piece polycarbonate lens and a separate nose cup to reduce internal lens fogging.

Features / Benefits

- Single lens for wide field of vision, allows for ready use of night vision aids
- Easy to don and doff
- Filter canister may be mounted on either cheek, or in the center of the mask
- Accepts Military C2A1 Nuclear, Biological, and/or Chemical (NBC) Filter canister
- Fully elastic six-point head harness

Specifications

One size

Weight

2 lbs. (approx.)

Accessories (available at extra cost)

- Spectacle kit
- Fit Test Kit with Bayonet-type Attachment Adapter
- Carry bag

Approvals

- NIOSH approved with the CP3N filter cartridge for CN riot control agent and P100 particulate efficiency
- Lens meets requirements of ANSI Z87.1 1989 (Face and Eye Protection)

Manufacturer / Distributor

3M Occupational Health and Safety Division St. Paul, MN 55133 (800) 896-4223 www.3m.com

SGE-



3M "SGE-400" Negative Pressure, Full Facepiece Respirator

Equipment Type:

Negative Pressure Full Facepiece Respirators

Make / Model:

MSA "Millennium Chemical-Biological Mask" and "Advantage 1000 CBA-RCA Gas Mask"

Description:

These are tight-fitting, negative pressure respirators, which utilize single or dual filter cartridges. Both masks feature a flexible, 1-piece polyurethane lens bonded to a Hycar rubber facepiece, fully elastic six-point head harness, separate nosecup to reduce internal lens fogging, and a mechanical speaking diaphragm.

Features / Benefits

- Lightweight
- Large single lens allows wide field of vision
- Portable, easy to don and doff

Specifications

• "Millennium Chemical and/or Biological (C/B) Mask" uses standard NATO-thread filter canister, and incorporates an integral drinking tube for fluid intake in contaminated environments

Weight

- Millennium 1.45 lbs. (approx.)
- Advantage 1000 1.25 lbs. (approx.)

Accessories (available at extra cost)

- Electrostatic Precipitator (ESP) Communications System: one-piece self-contained electronic speech projection device
- Polycarbonate Lens Outserts
- Spectacle kit
- Butyl-coated Nylon Hood (Millennium only)
- Black vinyl carrier

Approvals

- NIOSH approved with P-100 particulate efficiency rating for protection against CN riot-control agents
- "CBA-RCA" Filter Canister tested against nerve (GB), and cyanide (AC, CK) war agents. See accompanying chart

Manufacturer / Distributor

Mine Safety Appliances, Inc. Corporate Headquarters **Pittsburgh, PA 15230** (800) MSA-2222 www.MSAnet.com

Millennium



Advantage 1000



MSA "Millennium Chemical-Biological Mask" and "Advantage 1000 CBA-RCA Gas Mask" Negative Pressure, Full Facepiece Respirators

ESP Communications System



Performance of Millennium/Advantage 1000 CBA-RCA Canister

Performance of Millennium/Advantage 1000 CBA-RCA Canister Against Various Chemical Agents

MSA's Millennium/Advantage 1000 CBA-RCA Gas Mask Canister has been tested to determine how effectively it can protect an individual against a variety of chemical agents (in accordance with CASHPAC specified on page 3).

This canister contains a pleated high-efficiency (HEPA) filter to remove aerosols and solid particulates, and an impregnated activated carbon bed to adsorb (filter out) gases and liquid vapors. The HEPA filter removes nearly all of the tear gas agents from the inspired airstream. Any small amount of tear gas agent that passes through the filter is adsorbed by the impregnated activated carbon. Sarin (GB), DMMP (a Sarin simulant), and HCN are also effectively adsorbed by the carbon bed.

The challenge agents are listed in column 1; testing conditions and requirements in columns 2, 3, and 4; and test results in the last column. The Flow Rate for all challenge agents was 64 ±1 Lpm, except for CK, which was 32 Lpm. Temperature was 25 ±3 degrees Centigrade.

MSA Millennium/Advantage 1000 CBA-RCA Canister Testing Conditions and Results

| Challenge Agent | Challenge Concentration, mg/m² | Relative Humidity, % | Minimum Service Time Requirement | RESULTS: Service Time until Break, in minutes |
|----------------------------------------------------------------------------------------------------------------------------------|--------------------------------------|-------------------------|----------------------------------------|-----------------------------------------------------|
| Dimethyl Methylphosphonate (DMMP) (Sarin Simulant) Nerve agent | 1000 ±70 | 25 ±0.5 | 60 minutes until a 0.04 mg/m² break | >63 |
| Sarin (GB) Nerve agent | 1000 ±70 | 25 ±0.5 | 60 minutes until a 0.04 mg/m² break | >70 |
| Hydrogen Cyanide (AC) Blood agent | 550 ±25 | 50 ±3 | 30 minutes until a 5 mg/m³ break | >40 |
| ICN after Equilibration for 6 hours @ 25% RH | 550 ±25 | 50 ±3 | 30 minutes until a 5 mg/m³ break | >40 |
| HCN after Equilibration for 6 hours @ 85% RH | 550 ±25 | 50 ±3 | 30 minutes until a 5 mg/m³ break | >40 |
| ortho-Chlorobenzylidene Malononitrile (CS) Tear Gas Irritant | 23 ±8 | 50 ±3 | 480 minutes until a 0.4 mg/m³ break | >480 |
| CS after Equilibration for 6 hours @ 25% RH | 23 ±8 | 25 ±3 | 480 minutes until a 0.4 mg/m² break | >480 |
| CS after Equilibration for 6 hours @ 85% RH | 23 ±8 | 85 ±3 | 480 minutes until a 0.4 mg/m² break | >480 |
| o. Chloroacetophenone (CN) Tear Gas Irritant | 101 ±6 | 50 ±3 | 480 minutes until a 0.3 mg/m³ break | >480 |
| CN after Equilibration for 6 hours @ 25% RN | 101 ±6 | 25 ±3 | 480 minutes until a 0.3 mg/m³ break | >480 |
| CN after Equilibration for 6 hours @ 85% RH | 101 ±6 | 85 ±3 | 480 minutes until a 0.3 mg/m³ break | >480 |
| Oleoresin capsicum (C ₁₈ H ₂₇ NO ₃), the "active" ingredient of the OC tear gas Irritant | 0.036 | 50 ±3 | None | >480 until a 0.001 mg/m³ break |
| Cyanogen Chloride (CK) Blood agent | Approx. 4000 | 80 | None | 5 to 9 until an 8 mg/m³ break |

Note 1. CN testing was performed occurring to Mil C-51560 by Edgewood Assenal, Note 2. GB testing was performed by TNO Netberlands.

Equipment Type:

Negative Pressure, Full Facepiece Respirators

Make / Model:

MSA "Phalanx Alpha RCA Gas Mask"

Description:

This is a tight-fitting, negative pressure respirator utilizing a single filter cartridge that may be mounted on either cheek. The mask comprises a single polycarbonate lens mounted in a EPDM-rubber facepiece, separate nosecup to reduce internal lens fogging, integrated mechanical speaking diaphragm, and a fully elastic five-point head harness.

Features / Benefits

- Large, coated lens resists scratching, allows wide field of vision
- Easy to don and doff
- Ease of disassembly for cleaning and parts replacement
- Compact size for use with visored helmets

Specifications

N/A

Weight

1.6 lbs. (approx.)

Accessories (available at extra cost)

- ESP Communications System: one-piece self-contained electronic speech projection device
- Polycarbonate Lens Outserts
- Spectacle kit

Approvals

- NIOSH approved with P-100 particulate efficiency rating for protection against CN riot-control agents
- "CBA-RCA" Filter Canister tested against nerve (GB), and cyanide (AC, CK) war agents, see accompanying chart

Manufacturer / Distributor

Mine Safety Appliances, Inc. Corporate Headquarters Pittsburgh, PA 15230 (800) MSA-2222 www.MSAnet.com

Phalan



MSA "Phalanx Alpha RCA Gas Mask" Negative Pressure, Full Facepiece Respirator

ESP Communications System



Equipment Type:

Negative Pressure, Full Facepiece Respirators

Make / Model:

GIAT-CenTech Group, Respiratory Protection Kit (RPK)-1

Description:

RPK-1 is comprised of eight "ARF-C" full face, negative pressure respirators, 16 ABK/NBC filter cartridges, and eight each of the following: Carrying bags; cleaning sets; optical frames; carrying harnesses. The ARF-C is a tight-fitting, negative pressure respirator comprised of a single polycarbonate lens mounted in a polyurethane facepiece, single, centermount filter cartridge, separate nosecup to reduce internal lens fogging, integrated mechanical speaking diaphragm, and a fully elastic six-point head harness.

Features / Benefits

- Modular Kit approach simplifies acquisition, storage, accountability, and deployment
- Large lens provides wide field of vision
- Easy to don and doff

Specifications

N/A

Weight

1.25 lbs. (approx.)

Accessories (available at extra cost)

• External, one-piece self-contained electronic speech projection device

Approvals

- European EN-316 Standard
- NATO Standard

Manufacturer / Distributor

GIAT Industries / CenTech Group Arlington, VA 22203 (703) 525-4444

ARF-C





RPK-1

ARF-C Negative Pressure, Full Facepiece Respirator and GIAT-CenTech Group, Respiratory Protection Kit (RPK)-1

Equipment Type:

Powered Air Purifying Respirator (PAPR)

Make / Model:

3M, Breathe-Easy 7 (BE-7) Rapid Response Powered Air System (RRPAS)

Description:

This is a high filtration, powered, air purifying respirator. In the RRPAS configuration the system components include: Cordura-nylon carry bag, which converts to a vest supporting the back mounted, Breathe-Easy Turbo Unit Filter Pack with three filter canisters; a rechargeable Ni-Cd battery; and the BE-7 neoprene rubber, tight fitting, full face mask with breathing tube and clamps.

Features / Benefits

- Continuous flow of filtered air reduces pulmonary stress
- Easy to don and doff
- Large front mounted pockets on the vest, supports storage of other gear
- Velcro" fastened panels permit wearer identification
- Tight fitting BE-7 face piece provides negative pressure protection in the event of battery failure

Specifications

- Battery Life/ Re-chargeable Ni-Cd up to 8 hours continuous use, 4-hour recharge time, up to 500 discharge / recharge cycles
- Non-rechargeable lithium battery 12-hours of continuous use, 10-year shelf life
- Positive-pressure airflow @ 4 cfm minimum
- Adjustable head harness

Weight

10 lbs. (approx.)

Accessories (available at extra cost)

- Battery charger Fully automatic, operates on 120v-60Hz, three sizes available
- Single unit, 5-outlet unit and 10-outlet unit
- Spectacle kit

Approvals

- Facemask meets requirements of ANSI Z87.1-1979 (Eye and Face Protection)
- NIOSH approved with the AP-3 and AEP-3 Filter Cartridge

Manufacturer / Distributor

3M Occupational Health and Safety Division St. Paul, MN 55133 (800) 896-4223 www.3m.com



3M, Breathe-Easy 7 (BE-7) Rapid Response Powered Air System (RRPAS), Powered Air Purifying Respirator (PAPR)

Powered Air Purifying Respirator (PAPR)

Make / Model:

3M, Breathe-Easy 10 (BE-10) Positive-Pressure, Butyl Hood Configuration

Description:

This is a high filtration, powered, air purifying respirator. In the BE-10 configuration the system components include: The belt mounted Breathe-Easy Turbo Unit Filter Pack with three filter canisters, a rechargeable Ni-Cd or non-rechargeable Lithium battery, the BE-10 butyl rubber hood with elasticized neck damn with breathing tube and clamps.

Features / Benefits

- Continuous flow of filtered air reduces face-piece fogging and pulmonary stress
- Easy to don and doff
- Loose-fitting hood accommodates glasses and facial hair
- Butyl rubber Hood/Cape provides liquid splash protection

Specifications

- Battery Life / Rechargeable Ni-Cd up to 8 hours of continuous use, 4-hour recharge time, up to 500 discharge / recharge cycles
- Non-rechargeable lithium up to 12 hours of continuous use, 10-year shelf life
- Positive-pressure airflow @ 6 cfm
- Elasticized headband

Weight

8 lbs. (approx.)

Accessories (available at extra cost)

- Battery charger Fully automatic, operates on 120v-60Hz, three sizes available
- Single unit, 5-outlet unit and 10-outlet unit

Approvals

- Butyl rubber hood meets military standard MIL-C-51251A for resistance to chemical and biological warfare agents
- U.S. Army approved for Chemical Stockpile Emergency Preparedness Program (CSEPP) responders
- NIOSH / P100 particulate efficiency approvals with the AP-3 and AEP-3 Filter Cartridges

Manufacturer / Distributor

3M Occupational Health and Safety Division St. Paul, MN 55133 (800) 896-4223 www.3m.com



3M, Breath-Easy 10 (BE-10) Positive-Pressure, Butyl Hood Configuration, Powered Air Purifying Respirator (PAPR)

Powered Air Purifying Respirator (PAPR)

Make / Model:

MSA OptimAir 6A

Description:

This is a belt mounted, high filtration, powered air purifying respirator. The system components include a nylon mounting belt supporting the motor/blower module with rechargeable Ni-Cd battery, dual-OptiFilter cartridges, and Hycar rubber, Ultravue Facepiece with breathing tube.

Features / Benefits

- Continuous flow of filtered air reduces pulmonary stress, enhances wearer comfort
- Portable, easy to don and doff
- Tight fitting UltraView Facepiece provides negative pressure protection in the event of battery failure

Specifications

- Battery Life (Ni-Cd) 8 hours of continuous use, 6- to 16-hour recharge time, 250 discharge / recharge cycles
- Positive-pressure airflow at 4 cfm minimum
- Adjustable head harness

Weight

- 7.4 lbs. (approx.) with NiCd battery
- 5.4 lbs. (approx.) with Lithium battery

Accessories (available at extra cost)

- Battery charger Fully automatic, operates on 120v-60Hz, Single-Unit, or Single-Unit, Dual Rate
- Disposable lithium battery with 10-year shelf life and 12-hour service life
- Spectacle kit

Approvals

N/A

Manufacturer / Distributor

Mine Safety Appliances, Inc. Corporate Headquarters Pittsburgh, PA 15230 (800) MSA-2222 www.MSAnet.com



MSA OptimAir 6A, Powered Air Purifying Respirator (PAPR)

Powered Air Purifying Respirator (PAPR)

Make / Model:

GIAT – CenTech Group, Respiratory Protection Kits (RPK) 1 and 2

Description:

RPK-1 is comprised of eight "ARF-C" full face, negative pressure respirators, 16 ABK/NBC filter cartridges, and eight each of the following: carrying bags; cleaning sets; optical frames; harnesses. RPK-2 is comprised of 16 ABK/NBC filter cartridges, 4 dual-filter "CASU C420" single speed air distribution units, 4 air delivery hoses, 4 carrying belts and harnesses, and 8 lithium batteries. The PAPR system is comprised of the "ARF-C" respirator in combination with a "CASU C420" unit with breathing hose, carrying belt/ harness, lithium battery, and two ABK/NBC filter cartridges.

Features / Benefits

- Modular Kit approach simplifies acquisition, storage, accountability, and deployment
- Continuous flow of filtered air reduces pulmonary stress, enhances wearer comfort
- Portable, easy to don and doff
- Tight fitting facepiece provides negative pressure protection in the event of battery failure

Specifications

- Battery Life (Lithium) up to 12 hours of continuous use
- Positive-pressure airflow @ 115 lpm (4 cfm) minimum
- Adjustable head harness

Weight

Individual PAPR System – 10 lbs. (approx.)

Accessories (available at extra cost)

Contact manufacturer

Approvals

- NATO Standard xxxxxxxxx for ARF-C and CASU C420
- NIOSH approval for air flow and P3 particulate efficiency pending

Manufacturer / Distributor

GIAT Industries / CenTech Group Arlington, VA 22203 (703) 525-4444



GIAT – CenTech Group, Respiratory Protection Kits (RPKs) 1 and 2, Powered Air Purifying Respirator (PAPR)

Self-Contained Breathing Apparatus (SCBA)

Make / Model:

Bio-Marine/Neutronics, "Bio-Pack 240" Re-Breather

Description:

The Bio-Pak 240 SCBA is a semi-closed circuit system. The critical components are housed in a high impact Noryl Shell with harness assembly, which distributes the weight of the unit across the back and hips of the wearer. The breathing gas is circulated via a set of neoprene hoses and tight fitting, silicone rubber, positive-pressure face piece with head harness.

Features / Benefits

- Positive-pressure design combined with cooled breathing gas reduces pulmonary stress
- Metal threaded hose connections ensures durability

Specifications

- 21-cubic ft. capacity, composite oxygen cylinder rated @ 3,000 psi
- 4-hour duration @ constant flow rate of 1.78 liters per minute
- Brass tubing for high-pressure oxygen lines for added safety
- Audible alarm @ 20 to 25 percent remaining service life

Color

Gray

Size

N/A

Weight

39 lbs. fully charged with mask, hoses, and coolant

Shelf Life

All test and repair @ field / user level

Accessories

- Molded carrying case
- Spare oxygen bottles
- Lime Pak refill (44 gal. drum)
- Kevlar harness and straps available

Approvals

NIOSH / MSHA Certification TC-13F-185 and TC-13F-206 (use in IDLH atmospheres and low temperature / $+15\ ^{\circ}F$

Manufacturer / Distributor

BioMarine

Exton, PA 19341 (610) 524-8800

Bio-Pack





Bio-Marine/Neutronics, "Bio-Pack 240" Re-Breather, Self-Contained Breathing Apparatus (SCBA)

BIOPAK 240

A Breakthrough in Technology

Biomarine engineers have taken the closed-circuit principle to its most advanced state with the BioPak 240 four-hour positive-pressure self-contained breathing apparatus. The benefits of increased duration, low usage and maintenance cost, and greater protection create an unbeatable combination for mine rescue teams, hazardous materials teams and industrial workers who have to be in potentially dangerous atmospheres for extended periods of time.





BioPak 240 Development

Biomarine's resources were applied toward developing a product that addressed some of the shortcomings of other extended-duration SCBA. Mine rescue teams were canvassed for suggestions, and the result was a unit that meets the needs of the

Protection—The Real Reason for SCBA

While no breathing apparatus can absolutely guarantee that contaminants will not enter the mask, the addition of a slight pressure in the unit considerably raises the odds in the user's favor. By maintaining a positive pressure in the system, normal leakage should be outward, not inward. Protection factor testing has shown increases in safety of as much as 200 times the protection afforded by negative pressure equipment. The BioPak 240 can be used in atmospheres immediately dangerous to life or health (IDLH).

Why Choose a BioPak 240?

LONG DURATION

Regardless of the user's size or physical activity, the BioPak 240 delivers at least four hours duration.

EASE OF USE AND COMFORT

The BioPak 240 is simple to use. All controls and displays are easily accessible. The unit is extremely comfortable—the weight is carried on the user's hips instead of on the shoulders. The well-fitting silicone facepiece with built-in speaking diaphragm provides a durable, effective seal. Breathing resistance is low, reducing the wearer's fatigue.



PHOTO COURTESY OF ZINC CORP. OF AMERICA

COOLING SYSTEM

All closed-circuit SCBA generate or retain heat. The BioPak 240 incorporates a built-in "blue ice" coolant canister to keep the breathing gas temperature at a comfortable level. It's about a 30-second operation—just install the frozen canister into the unit



and it's ready to go. Other than normal cleaning, the coolant canister does not require maintenance; just refreeze for the next use.

LOW USAGE COST

The LimePak CO₂ absorbent is provided in bulk, reducing the per-use cost of the BioPak 240 to almost one-eighth the cost of other extended-duration units. There is no need for a separate training canister—one canister and one type of absorbent is used for everything.



CLEAR VISION

The BioPak 240 incorporates a Therma-Clear insert for clear, fog-free vision. Just apply the Final Solution anti-fog cloth to the insert, and forget about fogging for four hours.



· EASE OF MAINTENANCE

Unlike other extended-duration SCBA, the BioPak 240 requires virtually no preparation and set-up time. All repairs and testing can be done in the field, with a minimum of training. An optional videotape describes tests and procedures necessary to keep the BioPak 240 in perfect operating condition.



We Support You

If a disaster occurs (for example, a mine explosion), Biomarine maintains an emergency supply of units, parts and accessories to adequately supply the rescue effort. We're there when it counts.



PHOTO COURTESY OF ZINC CORP OF AMERI

H66-23

PERSONAL PROTECTIVE EQUIPMENT: NBC Protective Suits

Permeable Chemical Protective Clothing

Make / Model:

LANX Fabric Systems, Chemical Protective Overgarments and Undergarments

Description:

These are military-style, air-permeable, chemical protective garments that incorporate adsorption technology through the use of polymerically encapsulated, activated carbon to provide protection against chemical warfare agents. The proprietary fabric can be tailored into any clothing configuration, including gloves and boot-liners, in a variety of weights, to provide varying degrees of chemical protection. The fabric can be combined with a variety of shell fabrics, including nylon, cotton, NYCO-blends, and aramids, to provide liquid, vapor and aerosol protection as well as flame resistance.

Features / Benefits

- Lowbulk, air-permeable fabric reduces heat stress
- Stretchable fabric for better fit and enhanced wearer comfort
- Launderable

Specifications

- Fabric Weight 10 oz./yd.²
- **Fabric Thickness** 0.06 in.
- Breakthrough times for specified agents:
 - Mustard (H) Greater than 48 hours
 - Lewisite (L) Greater than 24 hours
 - Sarin (GB) Greater than 24 hours

Size

- S, M, L, and XL
- User specifications

Shelf Life

N/A

Accessories

None

Approvals

Type 1 fabric meets Military Specification MIL-U-44435, for Chemical Protective Undergarments (CPU)

Manufacturer / Distributor

LANX Fabric Systems **Newark, DE** (302) 451-3060



Overgarments



Undergarments

LANX Fabric Systems, Chemical Protective Overgarments and Undergarments, Permeable Chemical Protective Clothing

Permeable Chemical Protective Clothing

Make / Model:

GIAT-CenTech Group, "Tactical Operations Multipurpose Protective Suit" (TOMPS)

Description:

This is a military-style two-piece, front-closing, hooded chemical protective suit. The suit incorporates physical filtering, and adsorption technology through the inner layer of activated-carbon-impregnated foam to provide protection against liquid, vapor, and aerosol NBC agents. A complete ensemble requires the use of carbon-foam socks, appropriate footwear, butyl-rubber (or equivalent) outer gloves, and appropriate respiratory protection.

Features / Benefits

- Permeable, lightweight fabric reduces heat stress
- Cotton / Polyester blend outer fabric is water, oil resistant, and flame-resistant
- Available in modular kit form, for up to eight personnel per kit

Specifications

- **Fabric Weight** 190 gr/m²
- **Suit Weight** 5 lbs. (approx.)
- Breakthrough times for specified agents:
 - Mustard (H) Greater than 24 hours @ 10 g/m² challenge concentration

Size

Small, Medium, Large, and X-Large

Shelf Life

- 10 years in original packaging
- 30-days cumulative wear outside of contaminated area

Accessories

Modular Kits include utyl gloves and chemical protective socks

Approvals

- European Standard EN 468, Protective Clothing for use against liquid chemicals
- NATO approval for chemical warfare (C/W) protection, NSN 8415-99-001-9888/92

Manufacturer / Distributor

GIAT Industries / CenTech Group Arlington, VA 22203 (703) 525-4444

TOMPS Suit





GIAT-CenTech Group, "Tactical Operations Multipurpose Protective Suit" (TOMPS), Permeable Chemical Protective Clothing

Impermeable Chemical Protective Clothing

Make / Model:

Remploy / ADI Technologies, "TYVEK-F" Chemical Decontamination Suit

Description:

This is a military-style one-piece, front-closing, hooded chemical protective overgarment. The seams and fabric provide barrier protection against C/B warfare agents, and radioactive particulates. The barrier properties of the suit are not altered by exposure to standard chemical decontamination emulsions.

Features / Benefits

- Lightweight reduces risk of body overheating
- Non-woven, polymer coated fabric offers superior strength and abrasion resistance
- Low cost, disposable after use

Specifications

- Fabric Weight 110 gr/m²
- **Suit Weight** less than 1 lb.
- Permeation resistance, breakthrough times for specified agents:
 - Mustard (H) Greater than 48 hours
 - Lewisite (L) Greater than 24 hours
 - Sarin (GB) Greater than 24 hours

Size

Small, Medium, Large, and X-Large

Shelf Life

10 years in original vacuum pack

Accessories

None

Approvals

- European Standard EN 468, Protective Clothing for use against liquid chemicals
- NATO approval for C/W protection, NSN 8415-99-001-9888/92
- U.S. Army Technical Escort Unit (TEU) for emergency response to C/W agent incidents

Manufacturer / Distributor

Remploy United Kingdom ADI Technologies Inc. (U.S. Distributor) McLean, VA 22101 (703) 734-9626



Remploy / ADI Technologies, "TYVEK-F" Chemical Decontamination Suite, Impermeable Chemical Protective Clothing

Impermeable Chemical Protective Clothing

Make / Model:

Kappler CPF-3 Protective Overgarments

Description:

These are chemical protective suits, available in a wide range of styles. The seams and fabric provide barrier protection against C/B warfare agents, and radioactive particulates. The barrier properties of the suit are not altered by exposure to standard chemical decontamination emulsions. The CPF-3, one-piece suit with attached hood and sock boots is appropriate for use as the outer garment of a limited-use Level "B" or "C" ensemble.

Features / Benefits

- Multi-layer barrier film laminated to a polypropylene substrate offers superior strength and abrasion resistance
- Wide variety of styles available
- Choice of seam styles available
- Economical

Specifications

- **Fabric Weight** 4.6 oz. per square yd.
- Suit Weight 1.2 lbs. (approx. for one-piece hooded suit with sock boots)
- Permeation resistance / breakthrough times for specified agents:
 - Mustard (HD), Sarin (GB), Soman (GD), VX Greater than 2 hours

Shelf Life

5 years in original packaging

Color

Tan

Size

Small, Medium, Large, and X-Large

Accessories

Available as a component of the "Chem / Bio Response-Pak," available through TriCon Environmental, Inc. (see below)

Approvals

Chemical Warfare Agent permeation resistance testing in accordance with Military Standard (MIL-STD)-282, test methods 204.1.1, 206.1.2, and T-208

Manufacturer / Vendor

Kappler Protective Apparel and Fabrics Guntersville, AL (800) 633-2410 TriCon Environmental, Inc. Huntsville, AL (800) 854-4334

CPF-3 Protective Overgarments



Kappler CPF-3 Protective Overgarment, Impermeable Chemical Protective Clothing

Kappler Technical Services

1-800-633-2410

Chemical Warfare Agent Data vs. Chemical Protective Fabric (CPF) 3

| Chemical Warfare Agent | Breakthrough Time | Breakthrough Criteria | Test Method |
|--------------------------------------------------------------------------|----------------------|--------------------------|-------------|
| Bis (2-chloroethyl) sulfide (Mustard: HD) | > 2 hours | 4.0 ug/cm ² | 204.1.1 |
| Isopropyl methylfluorophosphonate (Sarin: GB) | > 2 hours | 1.25 ug/cm ² | 206.1.2 |
| Pinacolyl methylfluorophosphonate (Soman: GD) | > 2 hours | 0.5 ug/cm ² | T-208 |
| O-ethyl S-(2-diisopropylaminoethyl) methylphosphonothiolate (Nerve: VX) | > 2 hours | 0.5 ug/cm ² | T-208 |

Testing was conducted at Battelle Labs Inc. in accordance with MIL-STD-282, Filter Units. Protective Clothing, Gas-Mask Components, and Related Products: Performance-Test Methods.

WARNING: There are uses, environments, and chemicals for which these garments and/or fabrics are unsuitable. It is the responsibility of the user to review available data and verify that the garment and/or fabric is appropriate for the intended use and meets all specified Government and industry standards. **CAUTION:** Do not use for fire protection. Avoid open flame or intense heat.

Impermeable Chemical Protective Clothing

Make / Model:

Kappler NFPA 1991-Compliant, "Responder" Level "A" Suits (Limited Use)

Description:

These are one-piece, vapor-proof, fully encapsulating suits, available in front or rear entry style, with flat or expanded style back and include an aluminized fiberglass overcover. The suits incorporate covered exhaust valves, integral sock boots, and an internal waist belt. The NFPA 1991-compliant suits offer barrier protection against C/B warfare agents, and radioactive particulates combined with radiant heat protection.

Features / Benefits

Sewn and strapped seams, 48-in. gas-tight zipper with double overlap storm flaps, hook and loop closures, Polished 40-mil. PVC view window, three-layer glove set (Viton, Silvershield, and Kevlar knit)

Specifications / Physical Properties: (Kappler Disclaimer – "These are to be considered physical properties and **should not be construed** as specifications")

- FabricThickness: 21 mil.
- **Fabric Weight:** 8.5 oz. per square yd.
- Suit Weight: 12.5 lbs. (approx., with aluminized fiberglass overcover)
- Permeation resistance / breakthrough times for specified agents:
 - See accompanying chart

Shelf Life

Indefinite, contingent on periodic pressure test and visual inspection results

Color

Blue

Size

XS, S, M, L, and XL to 4XL

Accessories

Pressure Test Kit, Air Supply Pass-thru fittings

Approvals

- Chemical Warfare Agent permeation resistance testing in accordance with MIL-STD-282, test methods T-208 and T-209
- National Fire Protection Association (NFPA) 1991, 1994 Edition compliant

Manufacturer / Distributor

Kappler Protective Apparel and Fabrics Guntersville, AL 35976 (800) 633-2410

Responder Level A



Kappler NFPA 1991-Compliant, "Responder" Level "A" Suit (Limited Use), Impermeable Chemical Protective Clothing

Kappler Technical Services

1-800-633-2410

Chemical Warfare Agent Data vs. Responder Fabric and Double Taped Seams

| | Breakthrough | Breakthrough | |
|-------------------------------------|--------------|------------------------|-------------|
| Chemical Warfare Agent | Time | Criteria | Test Method |
| Bis (2-chloroethyl) sulfide | | | |
| | >480 minutes | 0.5 ug/cm^2 | T-209 |
| (Mustard: HD) | | | |
| Isopropyl methylfluorophosphonate | | | |
| | >480 minutes | 0.5 ug/cm^2 | T-208 |
| (Sarin: GB) | | | |
| Pinacolyl methylfluorophosphonate | | | |
| | >480 minutes | 0.5 ug/cm^2 | T-208 |
| (Soman: GD) | | | |
| Chlorovinylarsinedichloride | | 2 | |
| | >480 minutes | 0.5 ug/cm^2 | T-209 |
| (Lewisite: L) | | | |
| O-ethyl S-(2-diisopropylaminoethyl) | | | |
| methylphosphonothiolate | >480 minutes | 0.05 ng/cm^2 | T-208 |
| | | | |
| (Nerve: VX) | | | |

Testing was conducted at Battelle Labs Inc. in accordance with MIL-STD-282, Filter Units. Protective Clothing, Gas-Mask Components, and Related Products: Performance-Test Methods.

WARNING: There are uses, environments, and chemicals for which these garments and/or fabrics are unsuitable. It is the responsibility of the user to review available data and verify that the garment and/or fabric is appropriate for the intended use and meets all specified Government and industry standards. **CAUTION:** Do not use for fire protection. Avoid open flame or intense heat.

Impermeable Chemical Protective Clothing

Make / Model:

Kappler "Responder Plus" Level "A" Suit (Limited Use)

Description:

This is a one-piece, vapor-proof, fully encapsulating suit available in front or rear entry style, with flat or expanded style back. The suit incorporates covered exhaust valves, integral sock boots, and an internal waist belt. The suit fabric offers broad chemical hold-out with high abrasion resistance providing barrier protection against C/B warfare agents, and radioactive particulates.

Features / Benefits

- Soft-feel fabric enhances wearer comfort
- Sewn and strapped seams
- 48-in. gas-tight zipper with double overlap storm flaps, hook and loop closures
- Polished 40 mil. PVC view window
- Detachable Butyl gloves

Specifications / Physical Properties: (Kappler Disclaimer – "These are to be considered physical properties and **should not be construed** as specifications")

- Fabric Thickness: 29.8 mil.
- **Fabric Weight:** 9.6 oz. per square yd.
- **Suit Weight:** 8 lbs. (approx.)
- Permeation resistance / breakthrough times for specified agents:
 - See accompanying chart
- Shelf Life:
 - Indefinite, contingent on periodic pressure test and visual inspection results

Color

Hi-visibility orange

Size

XS, S, M, L, and XL to 4XL

Accessories

Pressure Test Kit, Air Supply Pass-Thru fittings

Approvals

- Chemical Warfare Agent permeation resistance testing in accordance with MIL-STD-282, test methods T-208 and T-209
- American Society for Testing and Materials (ASTM) F1001 Liquid and Gaseous Chemical Test Battery

Manufacturer / Distributor

Kappler Protective Apparel and Fabrics

Guntersville, AL 35976 (800) 633-2410

Impermeable Chemical Protective Clothing

Make / Model:

Kappler, "Responder Climate System Monitoring (CSM)" Level "A" Suit (Limited Use)

Description:

This is a one-piece, vapor-proof, fully encapsulating suit, available in front or rear entry style, with flat or expanded style back. The suit incorporates covered exhaust valves, integral sock boots, and an internal waist belt. The suit fabric offers broad chemical hold-out with high abrasion resistance providing barrier protection against C/B warfare agents, and radioactive particulates.

Features / Benefits

- 48-in. gas-tight zipper with double overlap storm flaps, hook, and loop closures
- Polished 40-mil. Polytron Verification Configuration (PVC) view window
- Sewn and strapped seams
- Detachable Butyl gloves

Specifications / Physical Properties: (Kappler Disclaimer – "These are to be considered physical properties and **should not be construed** as specifications")

- Fabric Thickness: 21 mil.
- **Fabric Weight:** 8.5 oz. per square yd.
- **Suit Weight:** 7 lbs. (approx.)
- Permeation resistance / breakthrough times for specified agents:
 - See accompanying chart

Color

Tan

Size

XS, S, M, L, and XL to 4XL

Shelf Life

Indefinite – contingent on periodic pressure test, and visual inspection results

Accessories

Pressure Test Kit; Air Supply Pass-Thru Fittings

Approvals

- Chemical Warfare Agent permeation resistance testing in accordance with MIL-STD-282, test methods T-208 and T-209
- ASTM F1001 Liquid and Gaseous Chemical Test Battery

Manufacturer

Kappler Protective Apparel and Fabrics Guntersville, AL 35976 (800) 633-2410

Responder



Kappler, "Responder CSM" Level "A" Suit (Limited Use), Impermeable Chemical Protective Clothing

Kappler Technical Services

1-800-633-2410

Chemical Warfare Agent Data vs. Responder Fabric and Double Taped Seams

| Chemical Warfare Agent | Breakthrough Time | Breakthrough Criteria | Test Method |
|-------------------------------------|----------------------|--------------------------|-------------|
| Bis (2-chloroethyl) sulfide | | | |
| | >480 minutes | 0.5 ug/cm^2 | T-209 |
| (Mustard: HD) | | | |
| Isopropyl methylfluorophosphonate | | | |
| | >480 minutes | 0.5 ug/cm^2 | T-208 |
| (Sarin: GB) | | | |
| Pinacolyl methylfluorophosphonate | | | |
| | >480 minutes | 0.5 ug/cm^2 | T-208 |
| (Soman: GD) | | | |
| Chlorovinylarsinedichloride | | | |
| | >480 minutes | 0.5 ug/cm^2 | T-209 |
| (Lewisite: L) | | | |
| O-ethyl S-(2-diisopropylaminoethyl) | | | |
| methylphosphonothiolate | >480 minutes | 0.05 ng/cm^2 | T-208 |
| | | | |
| (Nerve: VX) | | | |

Testing was conducted at Battelle Labs Inc. in accordance with MIL-STD-282, Filter Units. Protective Clothing, Gas-Mask Components, and Related Products: Performance-Test Methods.

WARNING: There are uses, environments, and chemicals for which these garments and/or fabrics are unsuitable. It is the responsibility of the user to review available data and verify that the garment and/or fabric is appropriate for the intended use and meets all specified Government and industry standards. **CAUTION:** Do not use for fire protection. Avoid open flame or intense heat.

Impermeable Chemical Protective Clothing (Limited Use)

Make / Model:

Manufactured from Du Pont "TYCHEM 9400" and "TYCHEM 10,000" fabric

Description:

TYCHEM 9400 and 10,000 fabric is non-woven, spunbonded olefin. These fabrics offer high abrasion resistance, combined with high levels of impermeability to CW agents, and provide barrier protection against biological warfare agents and radioactive particulates. CPC manufactured from these fabrics offer high levels of 'limited-use' protection, as a Level "A," "B," or "C" overgarment.

Features / Benefits

- Non-Woven, spun-bonded, polymer coated fabric offers superior strength and abrasion resistance
- High visibility color
- Economical (compared to reuseable suits)

Specifications

- Fabric Weight / Thickness / Color:
 - TYCHEM 9400 6.1 oz. per square yd. / 17 mils. / Yellow
 - TYCHEM 10,000 9.5 oz. per square yd. / 29 mils. / High Visibility, Lime Yellow
- Permeation resistance / breakthrough times for specified agents:
 - Fabric only
 - See accompanying chart

Size

N/A

Shelf Life

5 years

Accessories

N/A

Approvals

- Meets all fabric requirements under NFPA 1993
- U.S. MIL-STD-282, Method T-209, Protocols DN3 and DN4 for HD and L; Method T-208, Protocols DN5 and DN6 for GA, GB, GD, and VX (revised May 1997)

Manufacturer / Distributor

DuPont Corporation Wilmington, DE (800) 448-9835

Chemical Warfare Agent Permeating Testing DuPont Tychem® Chemical Barrier Fabric

For additional information on DuPont Tychem® fabrics, call 1-800-558-9329 (1-800-55-TYFAX).

13

Chemical Warfare Agent Permeation Testing **DuPont Tychem® Chemical Barrier Fabric**

Tychem® 10,000 Fabric

| Agent | Common Name | Protocol | Average Breakthrough Time | Minimum Detectable Permeation (µg/cm²) |
|-------|----------------|----------|------------------------------|-------------------------------------------|
| GA | Tabun | DN6 | >12 hrs. | <0.0002 |
| GB | Sarin | DN6 | >12 hrs. | <0.0002 |
| GD | Soman | DN6 | >12 hrs | <0.0002 |
| HD | Sulfur Mustard | DN4 | >12 hrs. | <0.2000 |
| L | Lewisite | DN4 | >12 hrs | <0.0120 |
| VX | | DN6 | >12 hrs. | <0.0002 |

Tychem® 9400 Fabric

| Agent | Common Name | Protocol | Average Breakthrough Time | Minimum Detectable Permeation (µg/cm²) |
|-------|----------------|------------|------------------------------|-------------------------------------------|
| GA | Tahun | DN5 | >12 hrs. | <0.0001 |
| GB | | DN5 DN6 | >12 hrs. >12 hrs. | <0.0001 <0.0001 |
| GD | Somen | DN5 | >12 hrs | <0.0001 |
| HD | Sulfur Mustard | DN3 DN4 | >12 hrs. >12 hrs | <0.1000 <0.1000 |
| L | Lowisite | DN3 DN4 | >12 fm 2 fms | <0.0060 <0.0042 |
| VX | | DN5 DN6 | >12 hrs >12 hrs | <0.0001 <0.0001 |

Fabric Test Protocols. All tests performed in triplicate for DuPont Nonwovens by an independent accredited laboratory.

Protocol DN3-MIL-STD-282, Method T-209 (HD) or modified for Lewisite, for 12 hours at 10 g/m².

Protocol DN4—MIL-STD-282, Method T-209 (HD) or modified for Lewisite, for 12 hours at 100 g/m² (total coverage).

Protocol DN5-MIL-STD-282, Method T-208 (GB) or modified for GA, GD and VX, for 12 hours at 10 g/m2.

Protocol DN6—MIL-STD-282, Method T-208 (GB) or modified for GA, GD and VX, for 12 hours at 100 g/m2 (total coverage).

Tychem* fabrics (Tyvek QC*, Tychem* SL, Tychem* 7500, Barricade*, Tychem* 9400 and Tychem* 10,000) should not be used around heat, flames, sparks or in potentially flammable or explosive environments.
 Garments made of Tychem* fabrics should have slip-resistant or antislip materials on the outer surface of boots, shoecovers or other garment surfaces in conditions where slipping could occur.

Impermeable Chemical Protective Clothing (TYCHEM 10,000 fabric)

Make / Model:

Lakeland Industries, NFPA 1991-Compliant, Level "A" Suits (Limited Use)

Description:

This is a one-piece, vapor-proof, fully encapsulating suit, available in front or rear entry style, with flat or expanded style back, and includes an aluminized overcover and Bata "HAZMAX" boots. The suit incorporates three covered exhaust valves, integral sock boots, and an internal waist belt. The NFPA 1991-compliant suits provide barrier protection against C/B warfare agents, and radioactive particulates combined with radiant heat protection.

Features / Benefits

- Choice of Aluminized PBI/Kevlar fabric blend, or Aluminized Fiberglass Over-cover
- Sealed interior and exterior seams
- Two-layer (10-mil. Teflon and 40-mil. polished PVC) face shield
- Three-layer glove set (Viton, Silver Shield, and Kevlar knit)

Specifications

- **Fabric Weight** 9.5 oz. per square yd.
- **Suit Weight** 25 lbs. (approx.)
- Permeation resistance / breakthrough times for specified agents:
 - Tychem 10,000 fabric only
 - See accompanying chart

Size

Small, Medium, Large, and X-Large to 5XL

Shelf Life

5 years in original packaging, if unused

Accessories

- Glove Quick-Disconnect Assembly; Air Supply Pass-Thru fittings; Pressure Test Kit
- Crotch and Elbow reinforcement available on request

Approvals

- Tychem 10,000 fabric meets all fabric requirements under NFPA 1993
- Tychem 10,000 fabric only U.S. MIL-STD-282, Protocols DN3 and DN4 for HD and L, and Protocols DN5 and DN6 for GA, GB, GD, and VX (revised May 1997) (Lakeland does not test finished suits under above protocols)
- NFPA 1991-certified for complete ensembles (certification is invalid if components are acquired separately)

Manufacturer / Distributor

Lakeland Industries

Decatur, AL (800) 645-9291



Lakeland Industries, NFPA 1991-Compliant, Level "A" Suit (Limited Use), Impermeable Chemical Protective Clothing (TYCHEM 10,000 fabric)

Impermeable Chemical Protective Clothing (TYCHEM 10,000 fabric)

Make / Model:

Lakeland Industries, "Deluxe" and "Standard" Level "A" Suits (Limited Use)

Description:

These are one-piece, vapor-proof, fully encapsulating suits, available in front or rear entry style, with flat or expanded style back. The suits incorporate two covered exhaust valves, integral sock boots and an internal waist belt. The suit fabric offers broad chemical holdout with high abrasion resistance providing barrier protection against C/B warfare agents, and radioactive particulates.

Features / Benefits

- Deluxe and Standard Suits 48-in. zipper with double storm flap and storage bag
- Deluxe Suit Two-layer (10-mil. Teflon and 40-mil. PVC) face shield; three-layer glove set (butyl, Silver Shield, and Kevlar knit); sealed interior and exterior seams
- Standard Suit Single-layer (40-mil. polished PVC) face shield, butyl gloves, and exterior seams sealed

Specifications

- **Fabric Weight** 9.5 oz. per square yd.
- **Suit Weight** 8 lbs. (approx.)
- Permeation resistance / breakthrough times for specified agents:
 - Tychem 10,000 fabric only
 - See accompanying chart

Size

Small, Medium, Large, and X-Large to 5XL

Shelf Life

5 years in original packaging, if unused

Accessories

- Glove Quick Disconnect Assembly; Air Supply Pass-Thru fittings; Pressure Test Kit
- Crotch and elbow reinforcement available on request
- Additional Level "A" and "B" suits and other styles available, contact manufacturer for complete list

Approvals

- Tychem 10,000 fabric meets all fabric requirements under NFPA 1993
- Tychem 10,000 fabric only U.S. MIL-STD-282, Protocols DN3 and DN4 for HD and L, and Protocols DN5 and DN6 for GA, GB, GD, and VX (revised May 1997) (Lakeland does not test finished suits under above protocols)

Manufacturer / Distributor

Lakeland Industries Decatur, AL (800) 645-9291









Lakeland Industries, "Deluxe" and "Standard" Level "A" Suits (Limited Use), Impermeable Chemical Protective Clothing (TYCHEM 10,000 Fabric)

Impermeable Chemical Protective Clothing (TYCHEM 9400 fabric)

Make / Model:

Lakeland Industries, Style 94400 and Style 94450 Level "B" Suits (Limited use)

Description:

These are one-piece, rear entry, encapsulating suits with elasticized wrists and attached sock boots with boot flaps. The Style 94400 is a flat back suit with one shrouded exhaust port; the Style 94450 is an expanded back suit with two shrouded exhaust ports. The suit fabric offers broad chemical holdout with high abrasion resistance providing barrier protection against C/B warfare agents, and radioactive particulates.

Features / Benefits

- 48-in. zipper with storm flap
- 20-mil. Polished PVC faceshield

Specifications

- **Fabric Weight** 6.1 oz. per square yd.
- **Suit Weight** 5 lbs. (approx.)
- Permeation resistance / breakthrough times for specified agents:
 - Tychem 9400 fabric only
 - See accompanying chart

Size

Small, Medium, Large, and X-Large to 5XL

Shelf Life

5 years in original packaging, if unused

Accessories

- Crotch and elbow reinforcement available on request.
- Additional Level "A" and "B" suits, and additional styles available, contact manufacturer for complete list

Approvals

- Tychem 9400 fabric meets all fabric requirements under NFPA 1993
- Tychem 9400 fabric only U.S. MIL-STD-282, Protocols DN3 and DN4 for HD and L, and Protocols DN5 and DN6 for GA, GB, GD, and VX (revised May 1997) (Lakeland does not test finished suits under above protocols)

Manufacturer / Distributor

Lakeland Industries Decatur, AL (800) 645-9291





Lakeland Industries, Style 94400 and Style 94450 Level "B" Suits (Limited Use), Impermeable Chemical Protective Clothing (TYCHEM 9400 Fabric)

Impermeable Chemical Protective Clothing (TYCHEM 10,000 fabric)

Make / Model:

MarMac NFPA 1991-Compliant "Commander Brigade" Level "A" Suit (Limited use)

Description:

This is a one-piece, vapor-proof, fully encapsulating suit, available in front or rear entry style with an expanded style back and includes an aluminized fiberglass overcover. The suit incorporates two covered exhaust valves, integral sock boots, and an internal waist belt. The NFPA 1991-compliant suit provides barrier protection against C/B warfare agents, and radioactive particulates combined with radiant heat protection.

Features / Benefits

- Sealed interior and exterior seams (interior suit)
- Wide view, four-layer face shield (inner suit 40-mil. PVC / 5-mil. Teflon / 20-mil. PVC, outer cover 20-mil. Lexan)
- Three-layer glove set (neoprene/ EVOH film/Kevlar knit)
- Knee wear pads (internal suit)

Specifications

- **Fabric Weight** 9.5 oz. per square yd.
- Suit Weight 13 lbs. (approx., includes inner suit and outer cover)
- Permeation resistance / breakthrough times for specified agents:
 - Tychem 10,000 fabric only
 - See accompanying chart

Size

Small, Medium, Large, X-Large, 2XL, and 3XL

Shelf Life

5 years in original packaging, if unused

Accessories

Glove Quick-Disconnect Assembly; Air Supply Pass-Thru fittings; Pressure Test Kit

Approvals

- Tychem 10,000 meets all fabric requirements under NFPA 1993
- Tychem 10,000 fabric only U.S. MIL-STD-282, Protocols DN3 and DN4 for HD and L, and Protocols DN5 and DN6 for GA, GB, GD, and VX (revised May 1997) (MarMac does not test finished suits under above protocols)
- NFPA 1991, 1994 edition compliant

Manufacturer / Distributor

MarMac Manufacturing, Inc McBee, SC (800) 845-6962



MarMac NFPA 1991-Compliant "Commander Brigade" Level "A" Suit (Limited Use), Impermeable Chemical Protective Clothing (TYCHEM 10,000 Fabric)

Impermeable Chemical Protective Clothing (TYCHEM 10,000 fabric)

Make / Model:

MarMac "Commander 10,000" Level "A" Suit (Limited use)

Description:

This is a one-piece, vapor-proof, fully encapsulating suit, available in front or rear entry style, with flat or expanded style back. The suit incorporates two covered exhaust valves, integral sock boots, and an internal waist belt. The suit fabric offers broad chemical holdout with high abrasion resistance providing barrier protection against C/B warfare agents, and radioactive particulates.

Features / Benefits

- Wide view, three-layer (40-mil. PVC / 5-mil. Teflon / 20-mil. PVC) Faceshield
- Double storm flap covered zipper
- Two-layer glove system (EVOH film + Butyl or Viton)
- Knee wear pads for added strength

Specifications

- **Fabric Weight** 9.5 oz. per square yd.
- Suit Weight -
- Permeation resistance / breakthrough times for specified agents:
 - Tyvek 10,000 fabric only
 - See accompanying chart

Size

Small, Medium, Large, X-Large, 2XL, and 3XL

Shelf Life

5 years in original packaging, if unused

Accessories

- Glove Quick-Disconnect Assembly; Air Supply Pass-Thru fittings; Pressure Test Kit
- "Commander" series also available in Tychem 9400 fabric, contact manufacturer for additional information

Approvals

- Tychem 10, 000 meets all fabric requirements under NFPA 1993
- Tychem 10,000 fabric only U.S. MIL-STD-282, Protocols DN3 and DN4 for HD and L, and Protocols DN5 and DN6 for GA, GB, GD, and VX (revised May 1997) (MarMac does not test finished suits under above protocols)

Manufacturer / Distributor

MarMac Manufacturing, Inc McBee, SC (800) 845-6962



MarMac "Commander 10,000" Level "A" Suit (Limited Use), Impermeable Chemical Protective Clothing (TYCHEM 10,000 Fabric)

Impermeable Chemical Protective Clothing (TYCHEM 10,000 and 9400 fabric)

Make / Model:

MarMac "Probe" Series, Level "B" Suits (Limited use)

Description:

These are one-piece, expanding back, encapsulating suits with two exhaust valves, elasticized wrists, and attached sock boots with boot flaps, available in front or rear entry style. Both fabric choices offer broad chemical holdout with high abrasion resistance providing barrier protection against C/B warfare agents, and radioactive particulates.

Features / Benefits

- Sealed external seams
- Wide view, three-layer (20-mil. PVC / 5-mil. Teflon / 5-mil. Polyester) face shield
- Double storm flap cover over PVC zip-lock zipper
- Knee wear pads for added strength

Specifications

- Fabric Weight 10,000 fabric / 9.5 oz. square yd. 9400 fabric / 6.1 oz. square yd.
- **Suit Weight** 8 lbs. (approx.)
- Permeation resistance / breakthrough times for specified agents:
 - Tyvek 10,000 fabric only
 - See accompanying chart

Size

Small, Medium, Large, X-Large, 2Xl, and 3Xl

Shelf Life

5 years in original packaging, if unused

Accessories

- Glove Ring Assemblies; Air Supply Pass-Thru fittings
- Additional suits and styles available, contact manufacturer for additional information

Approvals

- Tychem 10,000 and 9400 meets all fabric requirements under NFPA 1993
- Tychem 10,000 and 9400 fabrics only U.S. MIL-STD-282, Protocols DN3 and DN4 for HD and L, and Protocols DN5 and DN6 for GA, GB, GD, and VX (revised May 1997) (MarMac does not test finished suits under above protocols)

Manufacturer

MarMac Manufacturing, Inc. McBee, SC (800) 845-6962



MarMac "Probe" Series, Level "B" Suit (Limited Use), Impermeable Chemical Protective Clothing (TYCHEM 10,000 and 9400 Fabric)

Impermeable Chemical Protective Clothing (Re-Usable)

Make / Model:

MarMac "Commander UltraPro" Level "A" Suit

Description:

This is a one-piece, expanded back, vapor-proof, fully encapsulating suit, available in front or rear entry style. The suit incorporates two covered exhaust valves, integral sock boots, and an internal waist belt. The suit fabric offers broad chemical holdout with high abrasion resistance providing barrier protection against C/B warfare agents, and radioactive particulates.

Features / Benefits

- Flame resistant, single-layer, Teflon composite fabric
- Sealed interior and exterior seams
- Wide view multi-layer face shield (Inner 10-mil. FEP Teflon / Outer replaceable 20-mil. PVC overlay)
- Three-layer glove system (neoprene, EVOH film, and Kevlar knit)
- Knee wear pads for added strength
- Color-coded abrasion indicator layer

Specifications

- Fabric Weight –
- Suit Weight –
- Permeation resistance / breakthrough times for specified agents –

Size

Small, Medium, Large, X-Large, 2XL, and 3XL

Shelf Life

Indefinite – contingent on periodic visual inspection and pressure-test results

Accessories

Glove Quick-Disconnect Assembly; Air Supply Pass-Thru fittings; Pressure Test Kit

Approvals

- NFPA 1991, 1994 Edition certified
- Meets requirements of ASTM F1001

Manufacturer

MarMac Manufacturing, Inc. McBee, SC (800) 845-6962



MarMac "Commander UltraPro" Level "A" Suit, Impermeable Chemical Protective Clothing (Re-Usable)

Explosive Effects (Blast / Fragmentation) Protective Clothing

Make / Model:

Med-Eng Chemical / Biological SRS-5 Explosive Ordnance Disposal (EOD)-IEDD Protective Systems

Description:

The SRS-5 is a modular, full-body protective system comprised of trousers, jacket, groin protector, and helmet with choice of visor options. The SRS-5 system provides protection from explosives effects (thermal, blast overpressure, and fragmentation). The system is compatible with select C/B barrier suits (the manufacturer suggests the use of a permeable, carbon impregnated barrier garment and cooling suit), and respirators / self-contained breathing apparatus (SCBA).

Features / Benefits

Modular design combined with laminate construction offers inherent high level of protection against CW agent liquid and vapor

Specifications

N/A

Size

Small, Medium, Large, and X-Large

Weight

65 lbs. (approx.)

Shelf Life

N/A

Accessories

- Body Cooling Suit, BCS-3; Ballistic Add-ons; Hardwire and wireless Communications Sets; External Power supply and auxiliary batter-packs; Helmet-mounted search light
- Contact manufacturer for complete list of accessory / additional items

Approvals

Contact manufacturer for available test reports

Manufacturer / Distributor

Med-Eng Systems, Inc. **Ottawa, Ontario, Canada**(800) 644-9078 (North America only)



Med-Eng Chemical / Biological SRS-5 EOD-IEDD Protective Systems, Explosive Effects (Blast / Fragmentation) Protective Clothing

PERSONAL PROTECTIVE EQUIPMENT: PPE Accessories

Butyl Rubber Gloves

Description:

Butyl-rubber gloves exhibit the highest permeation resistance to gas or water vapors, the long cuff offers over-the-wrist protection. They remain flexible and sensitive at lower temperature. Curved fingers and hand design mean better fit and greater worker comfort. The gloves are durable and decontaminable.

Features / Benefits

- Variety of thickness' available
- Available with "Grip-Safe" textured palm

Specifications

N/A

Size

7, 8, 9, 10, and 11

Weight

N/A

Shelf Life

N/A

Accessories

None

Approvals

ASTM Standard F 739-91

Manufacturer / Distributor

North Safety Products Hand Protection Division Charleston, SC 29405 (803) 554-0660





Butyl Rubber Gloves

Silver Shield Gloves

Description:

Silver Shield gloves are highly impermeable to a wide range of solvents, acids, and bases. The long cuff design offers above-the-wrist protection. Made from "Norfoil" laminate they are lightweight, flexible, and readily conform to the hand.

Features / Benefits

- Low cost, disposable
- Ambidextrous, no pairing needed

Specifications

- Laminate construction
- 3-mil. thickness
- 14-in. length

Size

Small, Medium, and Large

Weight

N/A

Shelf Life

N/A

Accessories

None

Approvals

ASTM Standard F 739-91

Manufacturer / Distributor

North Safety Products Hand Protection Division Charleston, SC 29405 (803) 554-0660





HRR.E

Silver Shield Gloves

Nitrile Gloves

Description:

These are hypoallergenic, synthetic nitrile, form-fitting gloves. Designed to be worn as the first layer of gloves, they offer higher chemical, puncture, and abrasion resistance than rubber or latex gloves. They remain flexible at lower temperatures, which means increased wearer comfort.

Features / Benefits

- High level of dexterity with high tensile strength for superior resistance to cuts, snags and abrasion
- Ambidextrous
- Low cost
- Disposable

Specifications

- Available in two thickness':
 - "Standard" glove is 4 mils.
 - "Plus" glove is 8 mils.

Size:

- "Standard" available in M, Lg, and X-Lg
- "Plus" available in M and Lg

Weight

N/A

Shelf Life

N/A

Accessories

None

Approvals

Title 29 CFR 1910.1030, OSHA Standard for Occupational Exposure to Bloodborne Pathogens

Manufacturer / Distributor


Best Manufacturing Co. Menlo, GA 30731 (800) 241-0323 Associated Bag Co. Milwaukee, WI 53207 (800) 926-6100

Nitrile Gloves



H06-17

Nitrile Gloves



H56-18

Nitrile Gloves

Chemical Protective Boots

Make / Model:

"HAZMAX" Model 87011 and 87012

Description:

These are chemical protective knee-high boots suitable for wear with a level "A," "B," or "C" ensemble. These boots are manufactured from a PVC-polymer compound, which is 'injection molded,' resulting in a seamless boot offering full foot support, and protection from compression, impact, and underfoot penetration.

Features / Benefits

- Decontaminable
- Molded kickoff lug
- High visibility "Chemical Green" upper with yellow nonslip outer sole

Specifications

- 16-in. height
- Steel safety toe with full-length steel midsole

Size

Men's full sizes 6 through 15

Weight

Average 3 lbs. per pair

Shelf Life

N/A

Accessories

None

Approvals

- NFPA 1991 Chemical Permeation Resistance Requirements
- NFPA 1991 Flammability Resistance Test
- ANSI Z41 Class 75 for impact and compression
- ANSI Z41-PR for underfoot protection
- ANSI Z41-1983 for Dielectric Footwear

Manufacturer / Distributor

Bata Shoe Company, Inc. Belcamp, MD 20107 (410) 272-2000



"HAZMAX" Model 87011 and 87012, Chemical Protective Boots

Coveralls, Fire Resistant

Description:

These are one-piece, long sleeve, fire resistant coveralls suitable for wear as part of a Level "D" ensemble or as the inner layer of clothing for a Level "C" or "A" ensemble. They are zipper closed in the front using a two-way brass zipper.

Features / Benefits

- BI-swing back
- Two large breast and patch pockets
- Two pass through pockets with snap closure

Specifications

- NOMEX IIIA material.
- Colors Navy Blue, Royal Blue, Red, Orange, Gray, or Tan

Size

N/A

Weight

4.5 oz. per square yd. or 6 oz. per square yd.

Shelf Life

N/A

Accessories

None

Approvals

NFPA 1975 Standards on Work/Station uniforms for Firefighters, 1994 edition

Manufacturer / Distributor

Topps, Inc



Coveralls, Fire Resistant

Advent Safety Helmet

Description:

This safety helmet features a one-piece, multi-level shell. The helmet comfortably fits a large size range due to the Sure-Lock Ratchet system, 8-point suspension system, and adjustable chinstrap.

Features / Benefits

Multi-level construction, reduces shock to the neck and spine through front, back, top, and side impact
protection

Specifications

- Polycarbonate alloy outer shell
- Expanded polystyrene (EPS) foam core
- ABS inner shell
- Expanded polyethylene (EPE) crown

Colors

Blue, Red, Orange, White, or Yellow

Size

Adjustable

Weight

N/A

Shelf Life

N/A

Accessories

None

Approvals

- ANSI Z89.1-1986 Class A and B
- Cognizant Security Agency (CSA) Z94.1-1992 Class G

Manufacturer / Distributor

E.D. Bullard Company Cynthiana, KY

Advent Safety Helmet



H66-44

Advent Safety Helmet

Safety Glasses

Make / Model:

Encon Model 1920 Enfog

Description:

These are single lens, full frame safety glasses. They feature a one-piece wraparound frame, universal nose bridge, and antifog coating.

Features / Benefits

- High impact resistance
- Lightweight

Specifications

- Polycarbonate lens
- Two-post hinge mechanism incorporating a polycarbonate sideshield
- 99.9-percent UV filtration

Colors

- Frames / Clear
- Temple / Clear, White, Blue, Red, Purple, Yellow, and Orange
- Lens / Clear only

Size

N/A

Weight

N/A

Shelf Life

N/A

Accessories

None

Approvals

- ANSI Z87.1-1989
- CSA Z94.3

Manufacturer / Distributor

Encon Eye Protection Vallen Safety Supply Company Houston, TX 77040 (713) 462-8700

Encon Safety Glasses



H66-39

Encon Model 1920 Enfog Safety Glasses

Cooling Vest

Make / Model:

"Model 15"

Description:

This is a one-piece, six-pocket (three each front / back) vest. The cooling is achieved by inserting frozen "gel ice" strips in the pockets provided. It is designed to be slipped over the user's head and secured around the waist. It is intended to be worn under the Level "A" or "C" outer suit to help reduce the effects of heat stress.

Features / Benefits

- Low profile
- 4- to 6-hour duration

Specifications

- NOMEX or cotton construction
- "Thinsulate" brand insulation

Color

Blue

Size

Adult, teen, and child

Weight

8 to 12 lbs. fully loaded

Shelf Life

N/A

Accessories


Replacement gel strips

Approvals

N/A

Manufacturer / Distributor

ILC Dover, Inc. Frederica, DE 19946 (302) 335-3911



"Model 15" Cooling Vest

DETECTION EQUIPMENT

DETECTION EQUIPMENT:Detector Tubes and Tickets

Chemical Warfare Agent, Detection, and Identification (Reactive / PH Paper)

Make / Model:

Liquid Chemical Agent Detector (M8, M9, and 3-Way Paper)

Description:

M8, M9, and 3-way paper are used for the field detection and identification of liquid chemical warfare, nerve, and blister agent. M8 paper is issued in booklets of perforated, non-adhesive backed sheets; M9 and 3-way paper are issued in dispenser rolls of perforated, adhesive backed sheets. All are disposable, single use sheets.

Features / Benefits

- Easy to use
- Instantaneous response to indicate presence of liquid blister and nerve agent

Specifications

- M8 and 3-Way three-color response Yellow = G (non-persistent nerve agent), Dark Green to Black = VX (persistent nerve agent), Red = H (blister agent). Certain G-agents give red-brown color response that is intermediate between the typical H and G colors
- M9 light pink to red, single color response to indicate presence of liquid blister (H) and nerve agent (G and VX)

Color

N/A

Size

N/A

Weight

N/A

Shelf Life

- M8 and 3-way paper has indefinite shelf life, if kept dry
- M9 3 years from date of manufacture

Accessories

N/A

Approvals

N/A

Manufacturer / Distributor

Anachemia Canada Inc. Quebec, Canada (514) 489-5711 Truetech, Inc. Riverhead, NY 11901-2585 (516) 727-8600

M8/M9 Paper



H66-42

Liquid Chemical Agent Detector (M8, M9, and 3-Way Paper) – Chemical Warfare Agent, Detection, and Identification (Reactive / PH Paper)

Chemical Warfare Agent, Detection, and Identification, Colormetric Tubes / Tickets

Make / Model:

Chemical Agent Detector Kit, M18A2

Description:

This is a portable, manually operated detection kit used to detect and identify chemical warfare agent liquid and vapor in the field. It may be used for area or point sampling. The kit comprises a plastic carrying case, a set of laminated instruction cards, reagent solution/tablets, three plastic solution bottles (one green top, one blue top, one white top), enzyme tickets, Draeger-type tubes (yellow band, blue band, red band, green band, and white band), aspirator bulb, M8 paper, sample envelopes, and a pencil. There are no refill kits available to replace components of the kit expended in training or operations.

Features / Benefits

- Detects known chemical warfare agent liquids, vapors, and aerosols
- Easy to use

Specifications

- Reaction time less than 3 minutes per test, a complete series of tests may take up to 20 minutes
- Kit contains sufficient quantities for 25 complete series of tests
- Detection levels: GB 1.0 mg/m³, VX 0.1 mg/m³, H 0.5 mg/m³, HCN 8.0 mg/m³, L 10.0 mg/m³, CG 12.0 mg/m³, CX x.x mg/m³, AC/CK x.x mg/m³
- Provides for biological agent sampling; however, requires lab to analyze

Size

8 in. (L) x 6 in. (H) x 3 in. (W) (approx.)

Weight

2.5 lbs. (approx.)

Shelf Life

3 years from date of manufacture

Accessories

N/A

Approvals

Currently in use with U.S. military forces

Manufacturer / Distributor

Truetech, Inc. Riverhead, NY 11901-2585 (516) 727-8600

M-18 Kit





Chemical Agent Detector Kit, M18A2 –
Chemical Warfare Agent, Detection, and Identification, Colormetric Tubes / Tickets

Chemical Warfare Agent, Detection, and Identification, Colormetric Tubes / Tickets

Make / Model:

Chemical Agent Detector Kit, M256A1

Description:

This is a portable, manually operated kit used to detect and identify chemical warfare agent liquid and vapor. The kit consists of a plastic carrying case containing 12 individually wrapped and sealed sampler/detectors and a booklet of M8 paper. Unlike the M18 kit, the samplers in the M256 kit are ready to use and disposable after use; however, it does not allow the specificity of testing permitted by the M18 kit.

Features / Benefits

- Self-contained, requires no mixing of reagent or pre-testing of components
- Disposable after use

Specifications

- Detects the following nerve, blood, and blister agents shown, within 20 minutes:
 - G and V series nerve agents, H, HN, HD blister agents, L and CX blister agents, and AC and CK blood agents
- Kit contains 12 sampler detectors

Size

- 7 in. (L) x 3 in. (W) x 5 in. (H)
- Operating temperature range: -25 °F to 120 °F

Weight

1.2 lbs. (approx.)

Shelf Life

5 years from date of manufacture

Accessories

N/A

Approvals

Currently in use with U.S. military forces

Manufacturer

Acachemia Canada Inc. Quebec, Canada (514) 489-5711 Truetech, Inc. Riverhead, NY 11901-2585 (516) 727-8600

M256 Kit







 $Chemical\ Agent\ Detector\ Kit,\ M256A1-\\ Chemical\ Warfare\ Agent,\ Detection,\ and\ Identification,\ Colormetric\ Tubes\ /\ Tickets$

Chemical Warfare Agent, Detection, and Identification, Colormetric Tubes / Tickets

Make / Model:

MSA / AUER Chemical Warfare Agent Detector Tubes

Description:

These are "Draeger-type" colormetric tubes used to detect and identify chemical warfare, nerve, blister, and blood agent vapor.

Features / Benefits

- Does not require addition of reagent after sampling
- Rapid response time of xx minutes
- Easy to read "Go / No-Go" color change
- Disposable after use

Specifications

- Detects chemical agent vapors shown:
 - Part # 10007654 G-series and VX nerve agents
 - Part # 10007653 Mustard (H) blister agent
 - Part # 10007650 Lewisite (L) blister agent
 - Part # 10007652 Blister agent
 - Part # 10007651 Choking (CG/DP) and Blood (AC/CK) agent
- Specific detection threshold levels not available at time of printing

Size

Standard 7mm diameter tube

Weight

N/A

Shelf Life

N/A

Accessories

- MSA "Kwik-Draw" Deluxe Pump, "Gas-Tester II H" Pump, "Toximeter II" Automatic Pump
- Draeger "Model 31 Bellows Pump," or "Accuro Pump" are compatible with these tubes

Approvals

N/A

Manufacturer

Mine Safety Appliances Pittsburg, PA (800) MSA-2222 www.MSAnet.com



MSA/AUER Chemical Warfare Agent Detector Tubes – Chemical Warfare Agent, Detection, and Identification, Colormetric Tubes / Tickets

Chemical Warfare Agent, Detection, and Identification, Colormetric Tubes / Tickets

Make / Model:

Draeger "Civil Defense Set" (CDS), Chemical Warfare Agent Detection System

Description:

This is a manually operated, point or area air sampling system, comprising the following items: one carrying case; two Accuro hand pumps; five CDS Set I; five CDS Set V; one Tube opener; two Set openers; one Accuro pump instruction sheet; one spare parts kit; one 1-meter extension hose w/adapter; one reorder card; one laminated instruction card; and one indelible marker.

Features / Benefits

- Easy to use "Bellows" type pump, allows for one hand operation
- Corrosion resistant, needs no lubrication or special tools for maintenance
- Minimal user-training required
- Colormetric tubes are scaled to show concentration of gas or vapor present
- 5-minute response time

Specifications

- Uniform pump stroke delivers 100 cubic centimeters of sample air per stroke
- Each "CDS Set I" contains five, manifold-mounted tubes that detect the following agent / concentration:
 - AC (Hydrocyanic acid) @ 1 ppm; CG (Phosgene) @ 0.2 ppm; L (Lewisite) @ 3 mg/m³ (organic arsenic) 0.1 ppm (arsine); HN1, HN2, HN3 (Nitrogen Mustard) @ 1 mg/m³ (organic nitrogen compounds); H, HD (Sulfur Mustard) 1 mg/m³ (thioether)
- Each "CDS Set V" contains five, manifold-mounted tubes that detect the following agent / concentration:
 - Nerve Agents (phosphoric acid esters) @ 0.025 ppm; CG (Phosgene) @ 0.2 ppm; CK (Cyanogen Chloride) @ 0.25 ppm; Chlorine @ 0.2 ppm; HN1, HN2, HN3 (Sulfur Mustard) @ 1 mg/m³ (thioether)

Colors Weight

"Accuro" pump is blue with black handles N/A

Size Shelf Life

N/A Colormetric Tubes – dependent on type

Accessories

CDS-I and CDS-V sets may be acquired separately to replace expended sets

Approvals

N/A

Manufacturer / Distributor

National Draeger, Inc. Pittsburgh, PA 15230-0120 (412) 787-8383 Vallen Safety Supply Co. Houston, TX 77040 (713) 462-8700



Draeger "Civil Defense Set" (CDS), Chemical Warfare Agent Detection System – Chemical Warfare Agent, Detection, and Identification, Colormetric Tubes / Tickets

Chemical Warfare Agent, Detection, and Identification, Colormetric Tubes / Tickets

Make / Model:

Draeger Automatic Air Sampling / Gas (Vapor) Detection System

Description:

This is a battery powered, electro-mechanical pump housing that uses the Accuro hand pump. The system allows for automatic point or area air sampling. It is fully programmable and displays total number of strokes, number of strokes remaining and low battery condition. Other system components include the battery pack and charger, carrying strap, and Drager Colormetric Tubes.

Features / Benefits

- Eliminates manual operation of the "Accuro" hand pump
- Corrosion resistant, needs no lubrication or special tools for maintenance

Specifications

- LCD display
- 6-volt (rechargeable) battery

Colors

Blue and black

Size

7.5 in. x 10.75 in. x 4 in.

Weight

5.25 lbs.

Shelf Life

Tubes – dependent on type

Accessories

- Deluxe Kit includes the pump, tube opener, extension hose adapter, spare parts kit, screwdriver, 10 rubber caps, rigid carrying case, and permanent marker
- Soft Side Kit includes the pump, spare parts kit, screwdriver, permanent marker, and nylon carry pouch
- Other extension hose (3m, 10m, and 15m), hot air probe, and motor vehicle exhaust probe

Approvals

N/A

Manufacturer / Distributor

National Draeger, Inc. Pittsburgh, PA 15230-0120 (412) 787-8383 Vallen Safety Supply Co. Houston, TX 77040 (713) 462-8700

Accuro-2000



H66-28

Quantimeter



Draeger "Automatic Air Sampling/Gas (Vapor) Detection System Chemical Warfare Agent, Detection, and Identification, Colormetric Tubes / Tickets

DETECTION EQUIPMENT: IMS Instruments

Chemical Warfare Agent, Detection, and Identification, Ion Mobility Spectrometry

Make / Model:

Graseby, "Chemical Agent Monitor" (CAM)

Description:

The CAM is a military-style, hand-held point monitor designed to determine and indicate the hazard from chemical warfare nerve and blister agent vapor in the field. The system components include the detector unit, storage and carrying case, carrying harness, primary and spare lithium batteries, standoff collars, and simulant test kit.

Features / Benefits

- Back-lit Lightweight Chemical Detector (LCD) bar-graph display with audible alarm
- Rapid response and clearing times
- Two button operation
- Military Specification (Mil-Spec) design allows for decontamination with hot water / soap solutions

Specifications

- Air sample rate of 400 ml per minute
- Monitors for and discriminates between nerve (G–series and VX) and blister (HD and HN3) agent vapor at the minimum threshold shown:
 - GA, GB 0.1 mg/m³ in less than 30 seconds
 - $VX 0.4 \text{ mg/m}^3$ in less than 30 seconds
 - $HD 2.0 \text{ mg/m}^3 \text{ in less than } 15 \text{ seconds}$
- Battery life averages 10 hours @ 68 °F
- Operating temperature range / -13 °F to 113 °F

Size

15.5 in. (l) x 3.1 in. (W) x 5.6 in. (H)

Weight

3.9 lbs.

Shelf Life

- Requires Nuclear Regulatory Commission (NRC) license (or State equivalent) and life-cycle tracking for Nickle-63 source
- Requires annual maintenance and "swipe testing"

Accessories

N/A

Approvals

Currently in use with U.S. military forces

Manufacturer / Distributor

Graseby Dynamics, Limited Fairfax, VA 22030 (703) 218-0380



Graseby, "Chemical Agent Monitor" (CAM) Chemical Warfare Agent, Detection, and Identification, Ion Mobility Spectrometry

Chemical Warfare Agent, Detection, and Identification, Ion Mobility Spectrometry

Make / Model: Graseby, "Lightweight Chemical Detector" (LCD)

Description:

The LCD is a portable, hand-held area type air sampler that can detect subsymptomatic levels of select chemical warfare agents. It is user programmable in the field to detect additional live agents, or training/simulated agents.

Features / Benefits

- Easy-to-read LCD display
- Simultaneous detection of nerve, blister, blood, and choking agents
- Data-logging capability
- No NRC license requirement

Specifications

- Nonradioactive ionization source
- Detection of chemical agent vapor at the minimum threshold shown:
 - GD 0.1 mg/m³; GA 0.2 mg/m³; GB 0.5 mg/m³; HD 1.0 mg/m³
- 60-hours continuous operation with 4 AA alkaline batteries

Size

5.25 in. x 4.25 in. x 1.25 in. (approx.)

Weight

1.3 lbs. (with batteries)

Shelf Life

N/A

Accessories

N/A

Approvals

N/A

Manufacturer / Distributor

Graseby Dynamics, Limited Fairfax, VA 22030 (703) 218-0380



Graseby, "Lightweight Chemical Detector" (LCD) Chemical Warfare Agent, Detection, and Identification, Ion Mobility Spectrometry

Chemical Warfare Agent, Detection, and Identification, Ion Mobility Spectrometry Make / Model:

ETG "Advanced Portable Detector (APD) 2000," Chemical Agent Monitor / Detector

Description:

The APD-2000 is a hand-held, multi-purpose instrument designed to detect and discriminate between chemical warfare agent (nerve and blister) and riot control agent (Pepper Spray and Mace) in the field. System components include the detector unit; storage and carrying case; alkaline or rechargeable 'C' cell batteries, AC power adapter, and vehicle power adapter.

Features / Benefits

- Hand-held and lightweight
- User selectable "point source" or "area monitoring" operation
- Large, backlit LCD display, with audible alarm
- Rapid response and clear-down times
- Annual service available through the manufacturer

Specifications

- Detects and discriminates between nerve (G series and VX) and blister (HD, HN3, and L) agent vapor, and pepper spray and mace at the minimum threshold shown:
 - Nerve agent / 'G' series 15 ppb < 30 seconds, VX 4 ppb < 30 seconds
 - Blister agent / HD and HN3 300 ppb <15 seconds, L 200 ppb <15 seconds
- Operating temperature range: -22 °F to 126 °F
- Operates on internal batteries, external AC power, or 9 to 18 vdc

Size

4 in. x 3½ in. x 11 in.

Weight

6 lbs. (<3 kg) with batteries

Maintenance Requirements / Shelf Life

- Operator Service is equal to 5 minutes per 24 hours of operation
- Requires annual maintenance and "swipe testing" and NRC (or State equivalent) license for Nickle-63 source

Accessories

• Integrated radiation dosimeter; Data logging of all detection and monitoring events; Sample preservation; Remote telemetry via RF link

Approvals

N/A

Manufacturer / Distributor

Environmental Technologies Group, Inc. Baltimore, MD 21284-9840 (410) 321-5200

The threat could be as big as an all out

chemical attack...

...or as small as

pepper spray.

Either way, the APD 2000
(Advanced Portable Detector)
from Environmental Technologies
Group, Inc., offers full detection
and monitoring capabilities.
The first device of its kind to
detect chemical warfare agents,
pepper spray and mace, the
APD 2000 was designed for both
battlefield and civilian emergency
first response.



- Detects warfare agents, pepper spray and mace
- Selectable monitoring/detection
- Superior interferent resistance



Chem/Bio Marketing Manager 1400 Taylor Ave., P.O. Box 9840 Baltimore, MD 21284-9840 (410) 321-5370 Fax (410) 339-3123



ETG "Advanced Portable Detector (APD) – 2000," Chemical Agent Monitor / Detector Chemical Warfare Agent, Detection, and Identification, Ion Mobility Spectrometry

Chemical Warfare Agent, Detection, and Identification, Ion Mobility Spectrometry

Make / Model:

ETG "Fixed Site / Remote-Chemical Agent Detector" (FS/R-CAD)

Description:

The FS/R-CAD is a fixed or mobile sensor suite designed for use in the following applications: Force Protection, Perimeter Security, or Fixed-Site Monitoring as part of a remote detector network; mobile chemical warfare agent detection; and Hot Zone/Warm Zone monitoring. The FS/R-CAD employs the same detection and identification hardware and software as the APD-2000, with data transmission via RF-link, fiber-optic cable, or conventional hard-wire link.

Features / Benefits

- Custom enclosures available
- Simultaneous detection and identification of nerve and blister agents

Specifications

- Detects and discriminates between nerve (G –series and VX) and blister (HD, HN3 and L) agent vapor, and pepper spray and mace at the minimum threshold shown:
 - Nerve agent / 'G' series 15 ppb < 30 seconds, VX 4 ppb < 30 seconds
 - $-$ Blister agent / HD and HN3 $-\,300$ ppb <15 seconds, L-200 ppb <15 seconds
- Operating temperature range: -22 °F to 126 °F
- Operates on internal batteries, external AC power, or 9 to 18 vdc

Size

N/A

Weight

N/A

Shelf Life

Requires annual swipe testing and maintenance, and NRC (or State equivalent) license

Accessories

Radiation Detection capability, Meteorological sensors, Global Positioning System (GPS)

Approvals

N/A

Manufacturer / Distributor

Environmental Technologies Group, Inc. Baltimore, MD 2103 (410) 321-5370



ETG "Fixed Site / Remote-Chemical Agent Detector" (FS/R-CAD) – Chemical Warfare Agent, Detection, and Identification, Ion Mobility Spectrometry

DETECTION EQUIPMENT: Flame Spectrophotometry Instruments

Chemical Warfare Agent, Detection, and Identification, Flame Spectrophotometry

Make / Model:

AP2C, Automatic Chemical Warfare Agent Detector Kit

Description:

This is a portable, hand-held detector designed for field detection (area monitoring, or point sampling) of chemical warfare agent. The system comprises the AP2C Detector, S4PE Gun with SP4E Sample Scrapers (for liquid sampling), Hydrogen (H²) storage device, Sampling tube, Users Technical Manual, and Carry-bag/harness.

Features / Benefits

- Simple to use, one-button turn on
- Color coded instruction label on unit
- Rapid detection of agent
- Rapid clearing after detection of high level of agent
- Dual sensors provide simultaneous detection of nerve and blister agent vapor
- Easy-to-read Light Emitting Diode (LED) display

Specifications

- Detection of chemical agent vapor at the minimum threshold shown:
 - GB − 1.5 ppb, VX − 1.4 ppb, HD − 63 ppb
- 12 hours of continuous operation per battery / hydrogen container

Size

15.5 in. (L) x 3.1 in. (W) x 5.6 in. (H) (approx.)

Weight

4.4 lbs. with batteries

Shelf Life

N/A

Accessories

N/A

Approvals

N/A

Manufacturer / Distributor

GIAT-Centech Group Arlington, VA 22203 (703) 525-4444 Mine Safety Appliances Pittsburgh, PA 15230-0428 (800) 672-2222 www.MSAnet.com





AP2C, Automatic Chemical Warfare Agent Detector Kit – Chemical Warfare Agent, Detection, and Identification, Flame Spectrophotometry

DETECTION EQUIPMENT: SAW MiniCAD

Chemical Warfare Agent, Detection, and Identification, Sonic Acoustic Wave

Make / Model:

MicroSensor, Surface Acoustic Wave Miniature Chemical Agent Detector (SAW MiniCAD) Mk II

Description:

The SAW MiniCAD, Mk II is a portable, area type air sampler. It features automatic operation, low and high detection levels, audible alarm, and LED display for vital functions. Components include a "Pelican" storage container, users guide, external 6.3v sealed battery pack (with soft case), 120v, 60hz AC charger, detector unit with soft case and shoulder strap, RS232 serial communication cable (6-pin, male, mini DIN / DB-9 connector), and Vapor Confidence Check Source.

Features / Benefits

- Dual sensors provide simultaneous detection of nerve and blister agent vapor
- LED display for "System Check," "All Clear," "Low Battery," and agent alarm

Specifications

- Detection of chemical agent vapor at the minimum threshold shown:
 - $GD 0.1 \text{ mg/m}^3$, $GA 0.2 \text{ mg/m}^3$, $GB 0.5 \text{ mg/m}^3$, $HD 1.0 \text{ mg/m}^3$
- 2-minute warm up, 60-second response times
- Operates 3 hours on four internal lithium batteries, 6 hours on external battery
- Operating temp / humidity: +5 °C to +40 °C / 0% to 95% non-condensing

Size

5.25 in. x 4.25 in. x 1.25 in.

Weight

15 oz. with batteries

Shelf Life

5 years

Accessories


N/A

Approvals

N/A

Manufacturer / Distributor

Microsensor Systems, Inc. Bowling Green, KY 42103 (502) 745-0099



SAW MiniCAD Area Monitor

MicroSensor, SAW MiniCAD Mk II -Chemical Warfare Agent, Detection, and Identification, Sonic Acoustic Wave

DETECTION EQUIPMENT:Multi-Gas Monitors

Multi-Gas Monitoring

Make / Model:

BioSensors, "Cannonball 1700 Series" 4-Gas Detection Monitor

Description:

This is a portable monitor that can simultaneously monitor up to four atmospheric hazards. The monitor is housed in a ruggedly built, watertight body. The kit includes a storage / carrying case, shoulder strap, sampling draw system, rechargeable "gel-cell" battery and charger, and waterproof connector for an optional remote alarm. The monitor may also run on four alkaline D-cell batteries.

Features / Benefits

- Rapid response time
- Automatic backlit, LCD display
- One button control

Specifications

- User adjustable alarms, monitors
- Combustible gas when gas level exceeds preset alarm point
- Oxygen level High and Low levels
- Carbon monoxide and Hydrogen sulfide up to three alarm levels, ceiling, Time-Weighed Authority (TWA) and Short-Term Exposure Limit (STEL)
- Continuous sample draw pump with automatic liquid alarm and shutoff, prevents internal liquid contamination of meter
- Stores sample data in memory, downloadable to a PC

Size

91/4 in. x 61/2in. x 6 in.

Weight

6 lbs.

Shelf Life

N/A

Accessories (available at extra cost)

Calibration kits and gas, remote alarm, battery charger, and computer link kit

Approvals

N/A

Manufacturer / Distributor

Biosystems, Inc. Middlefield, CT 06455 (203) 344-1079 Vallen Safety Supply Co. Houston, TX 77040 (713) 462-8700

Multi-Gas Monitoring / Photo-Ionization Detector (PID)

Make / Model:

Rae Systems, "Multi-Rae"

Description:

This is a lightweight, compact, point or area sampler. The unit will simultaneously monitor oxygen level, combustible gas / Low Explosive Limit (LEL), presence of volatile organic compounds, and up to two user-defined toxic gases. Sample data may be stored internally and downloaded to a PC. The kit includes the Sampler / Meter with 10.6 ev lamp, wrist and shoulder strap, and belt loop and rechargeable Ni-Cd battery pack. The monitor may also run on four AA alkaline batteries.

Features / Benefits

- Programmable alarms
- Easy to operate, one button calibration
- Rugged construction

Specifications

- Internal sample pump, flow rate = 150 cc/min
- Monitors and displays: Combustible gas / 0 to 100% LEL, Oxygen level / 0 to 30
- VOCs 0 to 2000 ppm and up to two user defined toxic gases
- 90 dB audible and visual alarms
- One operation and two programming keys
- 10-hour continuous run time, may recharge and run simultaneously

Size

4.65 in. (L) x 3.0 in. (W) x 1.9 in. (H)

Weight

1 lb. with battery

Shelf Life

Lamp life varies with lamp voltage and usage, replacement intervals range from 2 to 24 months

Accessories

Calibration kits, External battery Charger, Automotive charging adapter, 15-ft. sample tube and remote sample probe, 11.6 eV lamp, Vibration Alarm, Earphone, Additional sensors, and Confined Space Kit

Approvals

N/A

Manufacturer / Distributor

Rae Systems Sunnyvale, CA 94086 (408) 481-4999 Vallen Safety Supply Co. Houston, TX 77040 (713) 462-8700

Photo-Ionization Detector (PID)



H66-7

Rae Systems, "Multi-Rae" – Multi-Gas Monitoring / Photo-Ionization Detector (PID)

DETECTION EQUIPMENT:Radiation Detection (Survey Instruments)

Radiation Detection / Survey

Make / Model:

Ludlum, Model 2241-2, Scaler / Ratemeter

Description:

The Model 2241-2 is a hand-held, general purpose survey meter. It is used to detect and measure gamma radiation and alpha-emitting contamination, and detect beta-emitting contamination. The system components include the Model 2241-2 Digital Scaler / Ratemeter, Model 44-6 Beta-Gamma G-M Detector, Model 44-9 Pancake G-M Detector, D-cell alkaline batteries, headphones, and a storage / carrying case.

Features / Benefits

- Aluminum housing with beige polyurethane enamel paint
- Backlit, four-digit LCD display includes "low battery" alert
- Two detector setup, allows interchangeable use of two probes with one meter
- Audible alert / alarm
- Automatic ranging, reduces potential for meter "flooding"

Specifications

- Accuracy is \pm 10% of true value with detector connected
- Operating temperature range is -4 °F to 122 °F
- Display in R/hr, Sv/h, cpm or cps when set in RATEMETER position
- Battery life / typically 200 hours with two "D" cell alkaline batteries

Size

Model 2241-2 – 6.5 in. (H) x 3.5 in. (W) x 8.5 in. (L)

Weight

3.5 lbs. with batteries

Shelf Life

N/A; however, periodic calibration with probes is required

Accessories

- Straight or coiled cables
- Headset
- Check source and holders
- Carrying case (Pelican or Airmold)
- Shoulder strap

Approvals

N/A

Manufacturer / Distributor

Ludlum Measurements, Inc. Sweetwater, TX 79556 (800) 622-0828 Atlantic Nuclear Corp. Canton, MA 02021 (617) 828-9118

Model 2241-2

DIGITAL SCALER/RATEMETER



H68-5

Ludlum, Model 2241-2, Scaler / Ratemeter – Radiation Detection / Survey

Radiation Detection / Survey

Make / Model:

Ludlum, Model 44-6 Geiger-Mueller (G-M) Detector Probe

Description:

The Model 44-6 G-M Detector is used for beta-gamma survey. It may be used for general area or point survey and personnel monitoring. It is used with the Model 2241-2 scaler / ratemeter. Operating voltage is supplied by the scaler / ratemeter. The probe is constructed of thin wall, polished stainless steel. The probe cable comes standard with a "C" series connector; other styles are available on request.

Features / Benefits

- Lightweight
- Durable construction and finish

Specifications (Model 44-6)

- G-M tube
- Beta-gamma survey sensitivity typically 1200 cpm/mR/hr
- Operating voltage is 900v

Size

6 in. (L) x 1.2 in. (D)

Weight

1 lb.

Shelf Life

N/A, however requires periodic calibration with scaler / ratemeter

Accessories

Straight or coiled cables

Approvals

N/A

Manufacturer / Distributor

Ludlum Measurements, Inc. Sweetwater, TX 79556 (800) 622-0828 Atlantic Nuclear Corp. Canton, MA 02021 (617) 828-9118

Model 44-6

THIN WALL G-M DETECTOR



Ludlum, Model 44-6 Geiger-Mueller (G-M) Detector Probe – Radiation Detection / Survey

Radiation Detection / Survey

Make / Model:

Ludlum, Model 44-9 Pancake G-M Detector Probe

Description:

The Model 44-9 Pancake G-M Detector is used primarily to detect and measure alpha contamination. It is used with the Model 2241-2 scaler / ratemeter. It may be used for general area or point survey, or for personnel monitoring. Operating voltage is supplied by the scaler / ratemeter. The pancake probe is constructed from painted aluminum.

Features / Benefits (Model 44-9)

Lightweight, easy to use

Specifications (Model 44-9)

- G-M tube
- Alpha, beta-gamma survey sensitivity typically 3300 cpm/mR/hr
- Operating voltage 900v

Size

1.8 in. (H) x 2.7 in. (W) x 10.7 in. (L)

Weight

1 lb.

Shelf Life

N/A, however requires periodic calibration with meter

Accessories

Straight or coiled cables

Approvals

N/A

Manufacturer / Distributor

Ludlum Measurements, Inc. Sweetwater, TX 79556 (800) 622-0828 Atlantic Nuclear Corp. Canton, MA 02021 (617) 828-9118

Model 44-9



Ludlum, Model 44-9 Pancake G-M Detector Probe – Radiation Detection / Survey

Radiation Detection / Survey

Make / Model:

MGP, DMC-2000 Series, Electronic Dosimeter

Description:

The DMC 2000 Series electronic dosimeter may be used as an area monitor, or to measure "whole body" exposure to ionizing radiation. It measures flat energy response to X-rays and gamma fields from 60 KeV to 6MeV, and linear response to dose rate fields of one mR/hr to 1000 mR/hr. The DMC 2000 series dosimeters may be used as a standalone device or integrated into a dosimetry system.

Features / Benefits

- Compact size, comparable to a personal pager
- 2,000 hour (1 year) continuous use with a standard calculator battery
- Audible and visual alarms
- Display units include uSv, mSv, mrem or cGy
- Water-proof, shock, drop, and vibration resistant

Specifications

- Accuracy is ± 5%
- Measurement range 0.1 mrem to 1,000 rem
- Operating temperature range is 14 °F to 122 °F
- Operating humidity up to 90%

Size

3.3 in. x 1.9 in. x 0.7 in.

Weight

2 oz. (w/battery)

Shelf Life

N/A

• Requires periodic calibration

Accessories / Related Items

- Dosimeters:
 - DMC 2000S X and Gamma (60keV to 6 MeV)
 - DMC 2000X X and Low-energy Gamma (20keV to 6MeV)
 - DMC 2000XB X, Low-energy Gamma, and Beta (Deep-dose/Shallow-dose, 20keV to 6MeV)
- Dosimeter Readers:
 - LDM 91 (Standalone Reader Station)
 - LDM 101 Dosimeter Reader (requires interface with PC)
 - LDM 2000 Dosimeter Reader (Pass-by data exchange)
- Dosimeter Calibrators:
 - CDM 21 Manual operation, or with optional robotic arm for automated operation
 - IRD 2000 Source check and recalibration

• Wireless Remote Monitors:

- WRMPlus Receiver Base station and teledosimetry applications
- WRM 2400 Receiver Base station and telemdosimetry applications
- CTM 2000 Transmitter Personal teledosimetry and area monitoring with DMC 2000
- CTM 2000E Transmitter Personal teledosimetry and area monitoring with DMC 2000

Area Monitor:

- AM 16 Remote area monitoring applications
- IM 2000 DMC 2000 interface module for remote monitoring applications
- EAM 100 DMC 2000 external alarm amplifier for remote area monitoring applications

• DOSISYS Software:

- DOSIMASS DMC 2000 and LDM 2000 maintenance and setup
- MPGAC Online access control
- DOSIVIEW Access control and geographic features/location tracking
- TELEVIEW Teledosimetry and area monitoring
- WINWRM Teledosimetry only
- DOSICAL Dosimeter calibration and reporting
- CDM 2000 Dosimeter calibration and reporting with the CDM 21

Approvals

N/A

Manufacturer / Distributor

MGP Instruments, Inc Smyrna, GA 30082 (770) 432-2744 www.mgpi.com

MGP Electronic Dosimetry System







MGP, DMC-2000 Series, Electronic Dosimeter – Radiation Detection / Survey

Radiation Detection / Survey

Make / Model:

Dosimeter Corporation, Self-Reading Dosimeter, Model # 862

Description:

The dosimeter is used to measure "whole body" exposure to ionizing radiation. It is a hermetically sealed, pocket type quartz-fiber electroscope. It features a cumulative, direct reading, 0-200 mRem scale. It requires the use of a charging device to charge the dosimeter and zero the scale prior to possible exposure.

Features / Benefits

- Easy to charge and use
- Accurate and rugged

Specifications

- Calibration Accuracy is ± 10% @ mid-scale
- Rate dependence is \pm 5% up to 10,000 R/h
- Electrical leakage = 2% per day maximum at room temperature
- Operating temperature range is -20 °F to 130 °F
- Operating humidity up to 98%

Size

½ in. (D) x 4½ in. (L)

Weight

1.25 oz.

Shelf Life

Requires periodic calibration and leak testing

Accessories

Dosimeter charger (Jordan Electronics Model V-750) with D-cell alkaline battery

Approvals


N/A

Manufacturer / Distributor

Dosimeter: Dosimeter Charger:

Dosimeter Corporation Atlantic Nuclear Corp.
Cincinnati, OH 45242 Canton, MA 02021
(513) 498-8100 (617) 828-9118

Self-Reading Dosimeter



H66-29

Dosimeter Corporation, Self-Reading Dosimeter, Model # 862 – Radiation Detection / Survey

Radiation Detection / Survey

Make / Model:

Landauer, Thermoluminescent Dosimeter (TLD), Type W8

Description:

This is a "wallet card" sized, film badge-type device used to measure "whole body" exposure to ionizing radiation. It is a completely passive device requiring no pre-use service by the user. The device contains two 100-percent lithium fluoride TLDs (chips), protected within a durable laminated finish. Requires lab or other specialized facility to "read" the badge and report exposure levels.

Features / Benefits

- Small, inconspicuous and unobtrusive, always on the job
- Clip to outer clothing or carry underneath

Specifications

Exposure range for standard badge is 30 m/rem to 1,000 rad for gamma and X radiation

Size

See attached

Weight

N/A

Shelf Life

Requires annual replacement of chip(s)

Accessories

- May be customized with up to four chips and indium foil
- May be customized to discriminate for beta radiation
- Custom logo's available


Approvals

NVLAP accredited for category 1 and 2 accidents (high and low energy photons)

Manufacturer / Distributor

Landauer, Inc. Glenwood, IL 60425-1586 (708) 755-7000

Thermoluminescent Dosimeter (TLD)



H86-32

Landauer, Thermoluminescent Dosimeter (TLD), Type W8 – Radiation Detection / Survey

DETECTION EQUIPMENT:Biological Agent – Detection and Identification

Biological Warfare Agent, Detection, and Identification

Make / Model:

Environmental Technologies Group, Inc. (ETG) Sensitive Membrane Antigen Reactive Test (SMART) Tickets

Description:

This kit provides capability for rapid detection of the generic presence of bacteria and biological toxins, and identification of specific bacteria and biological toxins. The kit consists of a sample collection component, a handheld instrument for detection of generic bacteria, a hand-held photometer for detection of proteinaceous toxins, and individual immuno-chemistry (SMART) tickets for specific biological agent, detection, and identification.

Features / Benefits

- Rapid response time: 5 to 15 minutes
- · Easy to use and read, minimal operator training required

Specifications

Identifies; Anthrax, Staphylococcal Enterotoxin 'B' (SEB), Y-Pestis (Plague), Botulinum Toxin (Bot Tox), Venezuelan Equine Encephalitis (VEE), and Brucella

Size

N/A

Weight

N/A

Shelf Life

N/A

Accessories

Training version available

Approvals

Similar technology currently in use with U.S. military forces

Manufacturer / Distributor

Environmental Technologies Group, Inc. Baltimore, MD 21284-9840 (410) 321-5200

For fast, precise biological detection,

ETG has just the ticket

Identify biological agents with convenient SMART™ Biological Warfare Agent Detection Tickets from Environmental Technologies Group, Inc. Each SMART ticket employs patented immuno-chemistry tests for specific biological agents.



SMART BIOLOGICAL WARFARE AGENT DETECTION TICKETS



Chem/Bio Marketing Manager 1400 Taylor Ave., P.O. Box 9840 Baltimore, MD 21284-9840 (410) 321-5370 Fax (410) 339-3123

- Quickly identifies agents
- Requires minimal training
- Prepare samples in minutes

ETG Sensitive Membrane Antigen Reactive Test (SMART) Tickets – Biological Warfare Agent, Detection, and Identification

Biological Warfare Agent, Detection, and Identification

Make / Model:

ETG "Bio-Detector"

Description:

The "Bio-Detector is a civilian variant of the Biological Integrated Detection System (BIDS), currently in use with the U.S. Armed Forces. The "Bio-Detector" provides an on-demand, inthe-field capability for the rapid detection and identification of specific biological warfare agents.

Features / Benefits

- Simultaneously detects up to eight different bio-agents
- 15-minute response time
- Audible and visual detection alarms
- Fully automated

Specifications

- 14-hour, continuous run time
- List of Bio-agents Identified available from the manufacturer
- Power requirements –

Size

24 in. x 24 in. x 30 in.

Weight

135 lbs. (approx.)

Shelf Life

N/A

Accessories


N/A

Approvals


Similar technology currently in use with U.S. military forces

Manufacturer / Distributor

Environmental Technologies Group, Inc. Baltimore, MD 21284-9840 (410) 321-5200



ETG "Bio-Detector" – Biological Warfare Agent, Detection, and Identification



ETG "Bio-Detector" – Biological Warfare Agent, Detection, and Identification (Continued)

DETECTION EQUIPMENT: NBC Event – Consequence Management

Weapons of Mass Destruction (WMD) / Nuclear, Biological, and/or Chemical (NBC) Agent, Detection / Survey

Make / Model:

NBC, Advanced Consequence Management System (ACMS)

Description:

The NBC-ACMS integrates NBC agent detection and identification equipment with hardware and software that performs personnel location tracking, asset/resource tracking, on-scene agency tracking, communications net information, vapor-plume modeling, and meteorological monitoring functions. The system provides the Incident Command System (ICS) / Incident Commander (IC) with the real-time, critical, information required for effective, WMD-NBC-event consequence management.

Features / Benefits

- Up to 2-kilometer (1.2-mile) range
- Broad agent coverage for NBC threats

Specifications

See appropriate page for each piece of equipment

Size

See appropriate page for each piece of equipment

Accessories / Related Items

- Minimum System Operating Requirements
- Optimal System Operating requirements
- Basic Components:
 - Hardware:
 - -- Two ea Advanced Portable Detector (APD)-2000 Chemical Warfare Agent Detectors
 - -- Five ea DMC 2000 Series electronic dosimeters
 - -- One ea LDM 2000 Dosimeter reader/personnel; marker
 - -- Four ea Wireless data extenders
 - Software:

Weight

See appropriate page for each piece of equipment

Shelf Life

See appropriate page for each piece of equipment


- Consequence Management Board (CMB) database
- -- "MIDAS Plume Generation" or "MIDAS AT Plume Generation, Tracking, and Modeling"
- Optional Hardware / Sensors
- CTM 2000 Transmitters
- WRM (*Plus* or 2400)
- SMART Tickets
- Automated Bio-Detector
- "MultiRad" Radiological Survey Package
- FPD 2000 Fixed Site Chemical Agent Detector
- Meteorological Sensors

Approvals

N/A

Manufacturer / Distributor

Environmental Technologies Group, Inc. Baltimore, MD 21284-9840 (410) 321-5370 MGP Instruments, Inc Smyrna, GA 30082 (770) 432-2744 Ocean Systems Engineering Corp. Oceanside, CA 92056 (760) 722-3222 www.mgpi.com



NBC, Advanced Consequence Management System (ACMS) – WMD/NBC Agent, Detection / Survey

DECONTAMINATION EQUIPMENT

DECONTAMINATION EQUIPMENT:Mass Decontamination Systems

Mass Casualty Decontamination Systems

Make / Model:

Advanced Containment Systems, Inc. (ACSI), Mass Casualty Decontamination System

Description:

Custom designed and built, transportable systems including trailer mounted systems ranging in size from 12 ft. to 60 ft., and standalone land-based or shipboard buildings / rooms. All systems offer custom electrical and plumbing designs, water and chemical resistant interiors, negative-pressure or positive-pressure ventilation, air and run-off filter systems, variety of holding tanks / bladders, and automatic injection systems incorporated with the shower systems for pre-mixed decontamination or soap solution.

Features / Benefits

- Self-contained system
- Rapidly deployable
- · Fixed-site or mobile, large-scale decontamination capability

Specifications

N/A

Size

Customer specifications from 12 ft. to 60 ft.

Weight

N/A

Shelf Life

N/A

Accessories


- Contact manufacturer for complete list of accessory items and manufacturing options
- Other design / application options include mobile command posts and mobile laboratories
- Modular Mass-casualty Decontamination System / components (air-beam shelters, Lp water heater, and vertical shower booms / nozzles)

Approvals

N/A

Manufacturer / Distributor

Advanced Containment Systems, Inc. **Houston, TX 77075** (800) 927-2271 www.acsi.com



ACSI, Mass Casualty Decontamination System

Mass Casualty Decontamination Systems

Make / Model:

MODEC, Mass Casualty Decontamination System

Description:

Standard or custom designed, mobile truck or trailer mounted platforms that can deliver up to 4,000 gallons per hour of heated decontamination solution and fresh rinse water. Various systems offer automated decontamination solution injector system, expandable foldouts with hose bibs and reels, interior and exterior decontamination stations, onboard fresh water storage and "gray-water" capture and holding tanks, and a contaminated clothing chute for packaging clothing.

Features / Benefits

- Self-contained system
- Rapidly deployable
- Large-scale decontamination capability

Specifications

Common system components include three-part chemical resistant, seamless, epoxy flooring; seamless, chemical resistant foam injected walls; aluminum gurney ramps; fold-down stainless steel tables; enclosed external modesty barriers; six-port manifold with cascade bottle air system; and climate controlled interior

Size

Various

Weight

Various

Shelf Life

N/A

Accessories

- Choice of external air-beam shelters; UV and HEPA H-VAC systems
- Contact manufacturer for list of current accessory items and available options

Approvals

N/A

Manufacturer / Distributor

MODEC, Inc. Denver, CO 80022 (303) 286-1546 www.massdecon.com



MODEC, Mass Casualty Decontamination System

DECONTAMINATION EQUIPMENT:

Modular Mass Decontamination Systems

Modular Mass Casualty Decontamination Systems

Make / Model:

TVI, Modular Mass Casualty Decontamination Systems

Description:

The TVI systems are lightweight, portable, and rapidly deployable. The system components may be stored and transported in a dedicated trailer, or prime mover (Truck, Squad, or Engine, etc.). The TVI system is comprised of the "Quik-Fire" diesel-fired flash hot water heater, user-selectable choice of PVC, copper, or stainless steel plumbing, high-pressure, quick disconnect hoses, hand operated Flex wands for supplemental decontamination, runoff containment berms, and a variety of user selectable, rigid articulating frame, rapid deployment shelters.

Features / Benefits

- Self-supporting, rapid setup shelters require no external power source, blowers, or air source
- All shelters offer side entry ports to reduce "wind-tunnel" effect
- All shelters may be "complexed" for custom configuration
- System is self-contained (except for water), for rapid deployment
- Quik-Fire heater includes 5 gal. onboard fuel, and quick-connect fuel bib for external fuel source

Specifications

- Quik-Fire Flash Heater: 30 cu. ft. (approx.); weighs 245 lbs. (approx.); heats 48 °F water to 86 °F at 950 gph flow; user adjustable fluid injection system preset for 0.5% bleach, may be reset for soap solutions or other solutions; requires #1 or #2 diesel fuel or kerosene.
- 4x4 Decontamination Shelter: 11.5 ft. interior width and various lengths (13 ft., 17 ft., 21 ft.) provides space for four simultaneous shower lines. Includes replaceable liner, interior curtains, fitted floor basin, and ports and hangers for shower hoses and booms.
- 2x6 Decontamination Shelter: 9.0 ft. interior width and various lengths (13 ft., 17 ft., 21 ft.) provides space for two simultaneous shower lines. Includes replaceable liner, interior curtains, fitted floor basin, and ports and hangers for shower hoses and booms.
- All TVI rigid frame shelters include a storage/transport bag, ground cloth, and anchor kit. All shelters provide high load rates and will remain standing with 30-percent frame damage (if properly anchored). Shelter fabric is fire retardant and UV stabilized for long life.

Weight

- Quik-Fire Flash Heater: 245 lbs. (approx.) on wheeled cart
- Shelters are dependent on width and length selected, average is 250 lbs. (approx.)

Shelf Life

N/A

Approvals

N/A


Manufacturer / Distributor

TVI Corporation Glendale, MD 80022

Accessories / Related Items

- High Efficiency Particulate Air (HEPA)rated air filtration system
- Positive-pressure or negative-pressure Heating, Ventilation, and Air-Conditioning (HVAC) system
- Specialty trailers w/generators and light sets
- Replacement parts for all system components
- Quick-erect scene management markers and barricades

(800) 598-9712 / www.tvicorp.com





TVI, Modular Mass Casualty Decontamination Systems

163

Modular Mass Casualty Decontamination Systems

Make / Model:

Reeves / MODEC, Modular Mass Casualty Decontamination System (MCDS)

Description:

The Reeves / MODEC modular system is portable and offers rapid setup capability. The basic system components include the Modular MCDS-unit; a Reeves "Airshelta," inflatable air-beam shelter with removable flooring for "berm-style" runoff containment; and two stainless steel booms with three shower heads on each boom.

Features / Benefits

- Self-contained system (except for water)
- Rapidly deployable
- Solid state electronic ignition in MCDS activates upon water flow
- Adjustable siphon jet injection system delivers bleach or soap decontamination solutions at constant dilution ratio
- Power vented exhaust vents towards the hot zone

Specifications

- MCDS configured to operate on LP gas, may be configured to operate on other fuels (diesel fuel, natural gas, or kerosene)
- Input pressure / fitting size up to 300 psi., 1½ fire hose fitting
- Output pressure 70 psi
- Output flow / temperature 530 gph @ 96 °F, adjustable up to 105 °F

Size

- MCDS 48 in. x 24 in. x 16 in.
- "Airshelta" 12 ft. x 12 in. x 9 ft.

Weight

- MCDS information not available at time of printing
- "Airshelta" 240 lbs. (approx.)

Shelf Life

N/A

Accessories

- Flow-thru Decontamination Brush
- Hand-held sprayer
- Accordion style roller systems for the decontamination of litter patients
- "Airsheltas" are available in larger sizes, and may be constructed to customer specifications
- Contact manufacturer for complete list of accessory / optional items that are available

Approvals

N/A

Manufacturer / Distributor

MODEC, Inc. **Denver, CO 80022**

(303) 286-1546 / www.massdecon.com

Reeves Manufacturing, Inc Fredrick, MD 21701 (800) 328-5563 / www.reevesmfg.com





Reeves / MODEC, Modular Mass Casualty Decontamination System (MCDS)

Modular Mass Casualty Decontamination Systems

Make / Model:

Air-Lite, Decontamination Shower and Hot Water Systems

Description:

This is a portable, easy to set up decontamination shower designed for use with hazardous materials (HAZMAT) or other rapid response teams. The shower framework and plumbing is PVC pipe with quick-connect unions, the shower outer fabric is urethane coated nylon. The shower breaks down for easy storage and transport in a nylon duffel bag.

The Hot-Water system is LP gas fired, designed to provide a maximum temperature increase of 90 °F, to a maximum temperature of 140 °F. The water supply is provided via a manifold reducer that allows use with standard fire apparatus or garden hose style intake.

Features / Benefits

- Shower System:
 - Includes nine adjustable jet-spray nozzles and one hand-held spray nozzle
 - Garden hose style connection with on/off valve
 - Double-door system allows for complexing of multiple units
- Stainless Steel Hot-Water System operates two decontamination shelters (wash and rinse)

Specifications

- Hot-Water System:
 - 6-hour run time (with two 5-gal LP bottles)
 - One 11/2in. and one garden hose intake (150-psi maximum)
 - Dual garden hose discharge outlets
 - Adjustable siphon injection system
 - Adjustable temperature control
 - 3-gpm flow rate

Size

- Shower Information unavailable at time of printing
- Hot-Water System 26 in. (W) x 28 in. (D) x 53 in. (H)

Weight

- Shower Information unavailable at time of printing
- Hot-Water System 255 lbs. (approx.)

Manufacturer / Distributor

Air-Lite Inflatables, Inc. **Haines, OR 77040** (713) 462-8700

Shelf Life

N/A

Accessories (available at extra cost)

- Runoff Collection Pool
- Decontamination elevation grid
- Variety of air-beam inflatable shelters

Approvals

N/A



Air-Lite, Decontamination Shower and Hot Water Systems – Modular Mass Casualty Decontamination Systems

HAZMAT Decontamination Shower

Make / Model:

Model HDS 6012

Description:

This is a portable, easy to set up decontamination shower designed for use with HAZMAT response teams. It is constructed of 1½in. diameter PVC pipe. The shower is designed to be placed in the RMC HAZMAT pool, which is intended to prevent contaminated runoff. The water supply is provided via a manifold reducer, which allows use with standard fire apparatus. The shower breaks down and stows in a canvas duffel bag.

Features / Benefits

- Color coded for ease of assembly
- Reusable

Specifications

N/A

Size

- Collapsed 48 in. x 18 in. x 8 in.
- Setup –

Weight

N/A

Shelf Life

N/A

Accessories (available at extra cost)

Manifold system with 11/2 in. adapter, decontamination pool, decontamination elevation grid, and shower curtain

Approvals

Currently in use with Washington Metropolitan Area (WMA) Metropolitan Medical Strike Team (MMST)

Manufacturer / Distributor

RMC Medical Philadelphia, PA 19154 (215) 824-4100 Vallen Safety Supply Co. Houston, TX 77040 (713) 462-8700



HAZMAT Decontamination Shower



Model HDS 6012 – HAZMAT Decontamination Shower

DECONTAMINATION EQUIPMENT:

Secondary/Individual Decontamination Systems and Supplies

Secondary / Individual / Equipment Decontamination

Make / Model:

M295 Kit, Decontamination, Individual Equipment

Description:

This is a personal use kit designed to decontaminate chemical warfare agents. The kit contains four hermetically sealed packets that contain one single use mitt. The mitt consists of non-woven fiber pad saturated with decontamination powder (synthetic resin mixture, E-555). After use, the mitt is disposed of as hazardous waste.

Features / Benefits

- Easy to use, blot, and wipe to remove agent
- Effective for decontamination of nerve and blister agents

Specifications

- Operating Temperature: -25 °F to 160 °F
- Mitt consists of non-woven fiber pad saturated with synthetic resin mixture, E-555 consisting of:
 - Sorptive resin
 - Strong acid ion-exchange resin
 - Strong base ion-exchange resin
 - Total water content is 25% by weight

Size

9 in. (L) x 5.5 in. (W) x 2 in. (T)

Weight

8 oz.

Shelf Life

N/A

Accessories

N/A

Approvals

Currently in use with U.S. military forces and WMA MMST

Manufacturer / Distributor

TrueTech, Inc. Riverhead, NY 11901-2585 (516) 727-8600

Secondary / Individual / Equipment Decontamination

Make / Model:

M291 Skin Decontamination Kit

Description:

This is a personal use kit for the removal of chemical warfare agents from the skin. The kit contains six individual decontamination packets containing a non-woven fiber applicator pad filled with decon powder (synthetic resin mixture, E-555). After use, the pads are disposed of as hazardous waste.

Features / Benefits

- Easy to use, blot, and wipe to remove agent
- Effective for skin decontamination of nerve and blister agents

Specifications

- Operating Temperature: -50 °F to 120 °F
- Decontamination pads consist of non-woven fiber pad saturated with synthetic resin mixture, E-555 consisting of:
 - Sorptive resin
 - Strong acid ion-exchange resin
 - Strong base ion-exchange resin
- Total water content is 25% by weight

Size

4.4 in. (H) x 9.3 in. (W) x 0.7 in. (T)

Weight

1.6 oz.

Shelf Life

N/A

Accessories

N/A

Approvals

Currently in use with U.S. military forces

Manufacturer / Distributor

TrueTech, Inc. Riverhead, NY 11901-2585 (516) 727-8600



H66-41

M291 Skin Decontamination Kit – Secondary / Individual / Equipment Decontamination

Secondary / Individual / Equipment Decontamination

Make / Model:

Canadian Reactive Skin Decontaminating Lotion (RSDL)

Description:

This is a broad spectrum, personal use, lotion-base decontaminant that is nonirritating to the skin and eyes, designed to decontaminate chemical warfare agents and select biological agents. RSDL is compatible with most skin types and, once applied, may be washed off with water.

Features / Benefits

- May be used as a wound decontaminant
- Offers protection as a "barrier cream" prior to exposure to chemical warfare agent
- Effective for skin decontamination of nerve and blister agents
- Residue and decontaminated agent is nontoxic

Specifications

- Self-indicating reactivity to agent, turns from transparent yellow to clear as reactive decontamination agent is consumed
- Single use pouch with one impregnated (RSDL 45g) foam pad
- 500-ml dispensing bottle

Size

N/A

Weight

N/A

Shelf Life

N/A

Accessories

N/A

Approvals

N/A

Manufacturer / Distributor

Therapex / E-Z-EM Canada, Inc.

O'Dell Engineering Ltd. Ontario, Canada (519) 740-8620 MODEC, Inc. Denver, CO (303) 286-1546

Secondary Decontamination System Supplies

Make / Model:

Description:

The items listed are intended to supplement the primary decontamination system and to provide limited decontamination capability where needed. The items may be procured locally.

| "Solo" Backpack 5-gal. sprayer | 15 ea. |
|--------------------------------|--------|
| 100-ft. garden type hoses | 10 ea. |
| Fan sprayers for hoses | 10 ea. |
| Sponges | 24 ea. |
| 20-in. brushes | 24 ea. |
| Fire hose reducing manifold | 3 ea. |

Features / Benefits

- Low cost
- Disposable
- Easily replaced

Specifications

N/A

Size

N/A

Weight

N/A

Shelf Life

N/A

Accessories

N/A

Approvals

N/A

Manufacturer / Distributor

Various, available locally



Secondary Decontamination System



Secondary Decontamination System Supplies

| Equipment Type: Decontaminant | | |
|-------------------------------|--|--|
| Make / Model: | | |
| | | |

Description:

The decontaminants listed are effective against nerve agent (Sarin and VX) and mustard (H), and are recommended for use on skin and equipment. The decontaminant includes dry and liquid soap used with water and various forms of active chlorine (liquid and dry).

Features / Benefits

- Readily available
- Easy to store and prepare

Specifications

- Chlorine content at least 5% by volume when mixed with water
- Bactricdal or germicidal soap, as appropriate

Size

Various

Weight

Various

Shelf Life

Rotate stock in accordance with manufacturer's specifications

Accessories

- Containers to mix and apply decontamination solution
- Brushes / sponges

Approvals

N/A

Manufacturer / Distributor

Various





Decontaminants

| Type of Agent | Effects | Onset | First Aid | Skin Decon | Detection |
|---------------------------------|------------------------------------------------------------|--------------------------------------|--------------------------------|-----------------------------------|-------------------------------------------------|
| Nerve: GA, GB, GD, GF, VX | Vapor: miosis, rhinorrhea, and dyspnea | Seconds | MARK I (1 to 3), Diazepam | M291, M258A1, bleach, water | M8A1 alarm, M256A1, M8, M9 papers, CAM |
| | Liquid: sweating, vomiting Both: convulsions, | Minutes to hours | | | |
| | apnea | | | | |
| Vesicants: H, HD, L | Erythema, blisters, irritation of eyes, cough, and dyspnea | Hours (immediate pain after L) | None | Same as above | M256A1, M8, M9 papers, CAM |
| Cyanide: AC, CK | Loss of conscious- ness, convulsions, and apnea | Seconds | None (nitrite and thiosulfate) | None usually needed | M256A1 |
| Pulmonary: CG (PFIB, HC) | Dyspnea, coughing | Hours | None | None usually needed | None |
| Riot Control: CS, CN | Burning, stinging of eyes, nose, airways, and skin | Seconds | None | Water | None |

| Decontaminant | Decontaminate | Mixing Procedures | Remarks |
|-------------------------------------------------------|--------------------------------------------------------------------------------------------|---------------------------------------------|-----------------------------------------------------------------|
| Sodium Hypochlorite Solution (Household Bleach) | Blister, V Agent, BIO, and CN (Tear) | 50/50 mixture of bleach and water | Flush with clear water after use |
| Sodium Carbonate (Washing Soda) | Vomiting, G, and CN | 10-lb soda to 12-gal water, 10% solution | Possible sources; commercial, laundry, and chemical firms |
| Soap and Water (Hot) | G and V on skin Liquid mustard agent, blood agent, arsenical, blister, tear, and vomiting | Mix thoroughly and keep hot | |
| M258A1 or M291 Decon Kit | G, V, mustards, arsenical, and blister | None | Only decontaminant available for paper |

Table 1. Decontaminants Using Active Chlorine (Refs. 1, 3, 4, 5, and 8)

| Decontaminant | Composition | Application |
|------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------|
| Household Bleach (Refs. 1, 3, 5, and 8) | 2-6 wt percent NaOCI in water | Skin and equipment |
| High-Test Hypochlorite (HTH) High-Test Bleach (HTB) (Refs. 1, 4, and 5) | Ca(OCI)CI + Ca(OCI) ₂ as a solid powder or as a 7-percent aqueous slurry (about 70-percent available chlorine) | Equipment or terrain |
| Super Tropical Bleach (STB) (Refs. 1, 4, 5, and 8) | Ca(OCI) ₂ + CaO as a solid powder or as 7-, 13-, 40-, and 70-wt percent aqueous slurries (30- to 35-percent available chlorine) | Equipment or terrain |
| Dutch Powder (Ref. 3) | Ca(OCI) ₂ +MgO | Skin and equipment |
| Activated Solution of Hypochlorite (ASH) (Refs. 3 and 5) | 0.5-percent Ca(OCI) ₂ + 0.5-percent NaH ₂ PO ₄ buffer + 0.05-percent detergent in water | Skin and equipment |
| Self-Limiting Activated Solution of Hypochlorite (SLASH) (Refs. 3 and 5) | 0.5-percent Ca(OCI) ₂ + 1.0- percent Na ₃ C ₆ H ₅ O ₇ + 0.2-percent H ₃ C ₆ H ₅ O ₇ + 0.05 percent detergent in water | Skin and equipment |
| German C8 Emulsion (Ref. 3) | Emulsion of 76-percent water, 15 percent perchloroethylene, 1-percent anionic surfactant, and 8-percent HTH | Equipment or terrain |
| Multipurpose Chemical/Biological Decontaminant (MCBD) (Ref. 3) | Microemulsion of 60-percent water, 5-percent HTH, 7-percent perchloroethylene, and the remainder cationic surfactant, cosurfactant, and sodium borate buffer (PH≈10) | Equipment or terrain |
| Chloramine (1,3-dichloro-5,5-dimethylhydantoin). The active component in Decontaminating Agent, Non-Corrosive (DANC) (Refs. 4, 5, and 8) | 6.25-percent chloramine in 1,1,2,2-tetrachloroethane | Equipment |

DECONTAMINATION EQUIPMENT:Decontamination System Accessory Items

Overpack Drums

Make / Model:

Description:

These are Polyethylene storage drums. The twist-on lids have a closed cell polyethylene gasket and are easily secured without tools. They provide a safe and economical means to transport and temporarily store containers of hazardous liquid or other contaminated material. They accept all steel drums up to 55-gal. size.

Features / Benefits

Durable and decontaminable

Specifications

N/A

Size

31 in. x 40¾in.

Weight

N/A

Shelf Life

N/A

Accessories

Steel drums up to 55 gal.

Approvals

Meets standards of Department of Transportation (DOT) 49, Code of Federal Regulations (CFR) 1733

Manufacturer / Distributor

Vallen Safety Supply Company **Houston, TX 77040** (713) 462 8700

Overpack Drums



Overpack Drums

Bags, Various

Make / Model:

Description:

- The "Personal Belongings" bags are for use by team members and personnel processing through the decontamination station.
- The "Cadaver" bags are to assist with mass fatality management.
- The "Laundry" bags and color-coded "polyethylene liners" are used primarily at the decontamination station for the collection of contaminated clothing and waste material.

Features / Benefits

- Personal Belongings: zipper style with write-on label
- Cadaver: contains fluids, center or side zipper closure, comes with 3-part ID tag and two color coded stickers for proper labeling
- Poly-Liners: sealed bottom, safe for landfill disposal or incineration
- Laundry: resistant to most organic solvents, impervious to most gases

Specifications

- Personal Belongings Clear with white "write-on" block, 4-mil. polyethylene, 8 in. (W) x 10 in. (L)
- Cadaver White, 10-mil. coated vinyl with dielectrically sealed seams. Maximum weight capacity of 350 lbs., 36 in. x 90 in.
- Poly-liners Red, 4 mil., with imprinted message "Contaminated Disposables," low density, free of heavy metals, 20 in. (W) x 13 in. (D) x 39 in. (L)
- Laundry Clear, 1¹/₄ mil., static free, completely soluble in water at 140 °F, can be heat sealed or with built-in ties, 36 in. (W) x 39 in. (L)

Size

See above

Weight

See above

Shelf Life

Indefinite

Accessories

N/A

Approvals


N/A

Manufacturer / Distributor

Associated Bag Co. Milwaukee, WI 53207-7120 (800) 926-4610

Personal Belongings Bags





Water Soluble Laundry Bags



Poly Liners



Poly Liners



Bags (Various)

Redress garments

Make / Model:

Description:

These are intended for use at the "clean" side of the decontamination station. The garments are to replace contaminated clothing removed by personnel as they enter the station and should be donned as they process out of the station. The garment types include disposable coveralls and paper hospital gowns.

Features / Benefits

- Lightweight, easy to store
- Range of sizes available
- Disposable after use

Specifications:

Appropriate for patient redress at decontamination station

Size

Various

Weight

N/A

Shelf Life

Indefinite

Accessories

Slippers / suitable foot covers

Approvals

N/A

Manufacturer / Distributor

Vallen Safety Supply Company **Houston, TX 77040** (713) 462-8700

MEDICAL / PHARMACEUTICAL SUPPLIES

Victim Resuscitation / Oxygen Administration

Make / Model:

Lifesaving Systems, "EM100 Oxylator" and "Oxylator Resus

Description:

The "EM100 Oxylator" resuscitation and inhalation system is designed to deliver a consistent volume of 100% O_2 or regular breathing air at positive pressure for patient ventilation. The "Oxylator Resuscitation Kit" provides the required equipment to resuscitate / ventilate up to four patients using compressed O_2 or compressed air from a conventional Self-Contained Breathing Apparatus (SCBA).

Features / Benefits

- Provides "hands-free" ventilation and ability to ventilate patients in confined spaces
- Easily cleaned and disinfected
- The system offers a "Manual" or "Continuous" mode of operation

Specifications

- "EM100 Oxylator":
 - Adjustable flow rate of 1.0 to 40 l/mn
 - Audible and visual inspiratory phase alarm
 - Operating temperature: -30 °C to 60 °C
 - Maximum deadspace = 20 ml
- Kit components include:
 - One Kit Bag; four "EM100 Oxylator" systems; four Coiled braided hose (length unspecified) with female DISS connection; two 50 ft. hoses with female DISS connection; two 30-ft. hoses with female DISS connection; four SCBA Regulators with gauges and DISS ports; four Ohio Female / DISS male quick-connects; four Ohio male / DISS male quick-connects

Size

N/A

Weight

N/A

Shelf Life

N/A

Accessories

Contact manufacturer for up-to-date accessory and related items list


Approvals

N/A

Manufacturer / Distributor

Lifesaving Systems, Inc Marietta, GA 30062 (770) 977-9907





Lifesaving Systems, "EM100 Oxylator" and "Oxylator Resuscitation Kit" Victim Resuscitation / Oxygen Administration

Mass Casualty / Trauma Supplies

Make / Model:

Superior Disaster Systems, Emergency Response Modules

Description:

Superior Disaster Systems offers a modular approach to mass casualty response readiness. The product line includes a variety of trailers prestocked with the requisite supplies to support a general medical response, to support mass casualty decontamination operations, or to support mass fatality management operations for a predetermined number of patients / casualties.

Features / Benefits

- Modular approach simplifies acquisition, storage, management, and deployment
- Components of the various modules include portable shelters, generators, light sets, tables and chairs, litters, dressings and bandages, tape, splints, face masks, decontamination showers, runoff collection pools, body / cadaver bags, etc.

Specifications

- Mass Casualty Medical Module (two sizes available)
- Hazardous Materials (HAZMAT) Decontamination Module
- Combination Decontamination and Medical Module
- Field Morgue Module

Size

N/A

Weight

N/A

Shelf Life

N/A

Accessories

Contact the manufacturer for complete list of modules, accessories, and related items

Approvals

N/A

Manufacturer / Distributor

Superior Disaster Systems, Inc. Thousand Oaks, CA 21075-6460 (805) 492-0393







Superior Disaster Systems, Emergency Response Modules – Mass Casualty / Trauma Supplies



Superior Disaster Systems, Emergency Response Modules – Mass Casualty / Trauma Supplies



Superior Disaster Systems, Emergency Response Modules – Mass Casualty / Trauma Supplies

Stretchers / Rescue Litter

Make / Model:

North American Rescue, "Raven" Rescue Litter

Description:

The "Raven" rescue litter was designed for emergency response teams. The frame is built of 1 1/8-in. aluminum alloy tubing with an automatic locking hinge, and incorporates ergonomically designed thermoplastic / fiberglass handles. The woven, mono-filament, polypropylene cover allows for fluid drainage and is easily replaced or decontaminated.

Features / Benefits

- Attached patient restraints with "Quick-Release" buckles
- All components are replaceable

Specifications

- Flame retardant materials
- Chemical resistant finish
- Payload rated at 300 lbs.

Size

- Folded 40 in. (L) x 7 in. (W) x 4 in. (H)
- Extended 81 in. (L) x 22 in. (W) x 6 in. (H)

Weight

11 lbs.

Shelf Life

N/A

Accessories

N/A

Approvals

N/A

Manufacturer / Distributor

North American Rescue Products, Inc. Elkridge, MD 21075-6460 (410) 799-2877

Raven Stretchers



North American Rescue, "Raven" Rescue Litter

Stretchers / Rescue Litter

Make / Model:

North American Rescue, "StingRay" Emergency Response Litter

Description:

The "StingRay" is a nonrigid litter designed for emergency response to a mass casualty event. The litter is constructed from a proprietary monofilament / polypropylene fabric that allows for maximum drainage of fluids while providing a nonslip surface for patient contact.

Features / Benefits

- Attached patient restraints with "Quick-Release" buckles
- Decontaminable fabric construction
- Built-in slots for carrying poles

Specifications

- Flame retardant and chemical resistant materials
- Six carrying points
- Folded 12 in. x 12 in. x 2 in.
- Extended 78 in. x 24 in.

Weight

3.3 lbs. (approx.)

Shelf Life

N/A

Accessories

N/A

Approvals

N/A

Manufacturer / Distributor

North American Rescue Products, Inc. Elkridge, MD 21075-6460 (410) 799-2877





North American Rescue, "StingRay" Emergency Response Litter

Specialized Pharmaceuticals

Make / Model:

Meridian Medical Technologies, Nerve Agent Antidote Auto-Injector

Description:

Meridian Medical technologies currently manufactures the MK1 Nerve Agent Auto-Kit (NAAK) and the CANA (Convulsive Antidote, Nerve Agent) Kit for the U.S. Department of Defense (DoD). Product improvements include a dual chamber auto-injector that combines atropine and 2PAM in a single injector.

Features / Benefits

- Easy to use, administers premeasured dose of antidote
- Small and compact, easy to transport

Specifications

- Mk1 NAAK comprises:
 - One "ATROPEN" auto-injector containing 2 mg of Atropine Sulfate Equivalent in 0.7 ml
 - One "COMBOPEN" auto-injector containing 600 mg Pralidoxime Chloride in 2 ml
- CANA (Diazepam) auto-injector contains 10 mg Diazepam in 2 ml

Size

N/A

Weight

N/A

Shelf Life

- MK1 NAAK / 5 years from date of manufacture
- CANA / 2 years from date of manufacture

Accessories

Training versions, without needles or drug product, are available

Approvals

U.S. Food and Drug Administration (FDA) approved for chemical agent treatment

Manufacturer / Distributor

Meridian Medical Technologies, Inc Columbia, MD (410) 309-1477



Meridian Medical Technologies, Nerve Agent Antidote Auto-Injector – Specialized Pharmaceuticals

Specialized Pharmaceuticals

Make / Model:

Taylor Pharmaceuticals, "Cyanide Antidote Package"

Description:

Formerly known as the "Pasadena Cyanide Kit," the Taylor Pharmaceuticals, "Cyanide Antidote Package" contains the required drugs to provide one complete treatment regimen for one (adult) cyanide exposure patient, or to provide the initial treatment regimen for two (adult) cyanide exposure patients, with the appropriate single use syringes and related items.

Features / Benefits

All required drugs and associated items required for administration are included in the kit

Specifications

- Each kit contains:
 - 2 ampules Sodium Nitrite injection, USP 300 mg in 10 ml water for injection
 - 2 vials Sodium thiosulfate injection, USP 12.5g in 50 ml water for injection
 - 12 ampules Amyl Nitrite Inhalants, USP 5 mg (0.3 ml)
 - 1 sterile 10-ml plastic disposable syringe with 22-ga needle
 - 1 sterile 60-ml plastic disposable syringe
 - 1 sterile disposable 20-ga needle
 - 1 stomach tube
 - 1 non-sterile 60-ml syringe
 - 1 tourniquet
 - 1 set of instructions for treatment of cyanide poisoning

Size

12 in. x 12 in. x 2 in. (approx.)

Weight

N/A

Shelf Life

2 years from date of manufacture

Accessories

N/A

Approvals

N/A

Manufacturer / Distributor

Taylor Pharmaceuticals San Clemente, CA 92673 (800) 223-9851

Specialized Pharmaceuticals

Make / Model:

Hope Pharmaceuticals, "Bulk" Organophosphate (Nerve Agent) and Cyanide Antidote

Description:

Hope Pharmaceuticals manufactures a variety of drug products. Specific products of interest include bulk quantity packaging of Sodium Nitrite and Sodium Thiosulfate for the treatment of cyanide poisoning, and Atropine Sulfate in soluble tablet form for reconstitution and intramuscular injection to treat organophosphate (nerve agent) poisoning.

Features / Benefits

- "Bulk" purchase of specialized drugs may offer overall cost reduction
- Acquisition of "bulk" quantities of specialized drugs offers flexibility in response planning and management for a mass casualty event

Specifications

- Sodium Nitrite Injection, USP 30 mg/ml, available in 10-ml ampules, two ampules per package
- Sodium Thiosulfate Injection, USP 12.5g in 50-ml vials, one vial per package
- Atropine Sulfate in 0.4-mg tablet form for reconstitution

Size

N/A

Weight

N/A

Shelf Life

2 years from date of manufacture for the drugs identified above

Accessories

N/A

Approvals

N/A

Manufacturer / Distributor

Hope Pharmaceuticals Scottsdale, AZ (602) 607-1970

GENERAL SUPPORT

Jumbo Turnout Bag

Make / Model:

Model P-BG065

Description:

This is a large and roomy equipment bag. It is constructed of 1,000 denier nylon with wrap-around reinforcing webbing and two webbed handles. It features a two-way zipper closure for easy access.

Features / Benefits

- Water resistant
- 5-in. "Maltese Cross" emblem screen printed on front of bag

Colors

Red with black handles and emblem

Specifications

- Single large compartment
- Zip closed exterior pocket

Size

16.5 in. (H) x 29.25 in. (W) x 12.5 in. (D)

Weight

N/A

Shelf Life

N/A

Accessories

None

Approvals

Currently in use with Washington Metropolitan Area (WMA) Metropolitan Medical Strike Team (MMST)

Manufacturer / Distributor

Gall's Inc. Lexington, KY 40555-4308 (800) 477-4766

Gall's Jumbo Turnout Bag



H66-40

Model P-BG065, Jumbo Turnout Bag

Command / Equipment Trailers

Make / Model:

Description:

Equipment trailers allow Metropolitan Medical Response System (MMRS) equipment to be stored in a centralized location. The components include a fully enclosed and lockable body on an appropriate frame / axle combination with hitch, safety cables or chains, lighting package (interior / exterior), level indicators, and leveling system.

Features / Benefits

- Allows for secure storage of equipment
- Provides centralized issue point of equipment
- Allows for rapid response / deployment

Specifications

N/A

Color

Paint scheme to meet user requirements

Size

User selectable

Weight

N/A

Shelf Life

N/A

Accessories

Tandem axle, enclosed cargo body with side and rear lockable doors, appropriate lighting (interior / exterior), shelving or other storage, work station with 110v electrical receptacles, external power supply receptacle(s), and heating / cooling capability, as required

Approvals

- WMA MMST is using an 18-ft., 7,000-lbs. GVW "Haul Mark" trailer
- AMA MMST is using a 20-ft., 7,700-lbs. GVW "Wells Cargo" trailer

Manufacturer / Distributor

Local Purchase



20-ft. Equipment Trailer



Equipment Trailer

Command / Equipment Trailers

207

Field Support Packaging

Make / Model:

Description:

These are reusable shipping containers, 17-gal. size, Model # 1BZ16. They are used to store and protect various items of equipment used by the MMST. They can be nested when empty, stacked when full. These are low cost and decontaminable, if required.

Features / Benefits

- Lightweight
- Lockable
- Easy to color code for identification of contents

Specifications

- Injection molded polyethylene
- Two-piece interlocking lid
- Stainless Steel hinge pins

Color

Gray

Size

211/ain. (L) x 15 in. (W) x 17 in. (H), 3.17 cubic ft.

Weight

9 lbs. (approx.)

Shelf Life

N/A

Accessories

Variety of sizes available

Approvals

Currently in use by AMA MMST

Manufacturer / Distributor

Ackro-Mills

Grainger Industrial /
Commercial Equipment & Supplies
345 Locations Nationwide
(800) 323-0620

AKRO - Mills Storage Container



H66-30a

Field Support Packaging

Field Support Packaging

Make / Model:

"Pelican" Cases

Description:

These are rigid equipment cases made from structural resin with foam packing inserts. The cases incorporate a neoprene O-ring seal and purge valve making them watertight and airtight. They are used to store and transport sensitive items of MMRS equipment.

Features / Benefits

- Watertight and airtight to 30 fsw
- Locking flanges and multiple latches

Specifications

N/A

Size

Various, see chart

Weight

Various, see chart

Shelf Life

N/A

Accessories


- Custom colors and inserts
- Personalized nameplate engraving and private logs available

Approvals

N/A

Manufacturer / Distributor

Pelican Products, Inc. Torrance, CA 90505 (310) 326-4700



H66-10

"Pelican" Cases, Field Support Packaging

PELICAN PROTECTOR® EQUIPMENT CASES

Watertight and airtight to 30 feet for the ultimate in protection. Constructed of light-weight space age structural resin with a neoprene o-ring seal and exclusive purge valve. Supplied complete with pre-scored pick n° pluck foam or padded dividers (no cutting required). Also includes locking flanges, massive multiple latches for absolute security and a comfortable molded fold down handle. Integral Personalized Nameplate engraved at optional cost. Custom colors and foam inserts fabricated on special order. Private label logos available. #1700 and #1750 supplied with regular foam, pick & pluck is not available.





1100 GUARD BOX® (NO FOAM)

Outside Dimensions: 8½" L x 6½"W x 3½" D. Inside Dimensions: 7½" L x 4¾" W x 3½" D. Colors: Black, Blue, Orange, Yellow, Lime Green, Hot Pink. Shipping Weight: 1 lb. 2 oz. Master Pack 6 ea.

#1100 F GUARD BOX® (With Foam)





1150 CASE

Outside Dimensions: 91/4" L x 73/4" W x 41/2" D. Inside Dimensions: 8½" L x 6"W x 3¾" D. Fold down handle. Colors: Black, Blue, Silver, Yellow, Hot Pink, Lime Green. Shipping Weight: 2.5 lbs. Master Pack 12 ea.



1200 CASE™

103/4"L Outside Dimensions: 103/4" L x 93/4"W x 5" D. Inside Dimensions: 9¹/₂" L x 7¹/₂" W x 4¹/₂" D. Fold down handle. Colors: Black, Silver, Yellow. Shipping Weight: 3 lbs. **Master Pack 6 ea.**




1300 CASE™

Outside Dimensions: 103/4" L x 93/4"W x 7" D. Inside Dimensio 91/2" L x 71/2" W x 61/2" D. Fold down handle. Colors: Black. Silver. Yellow. Shipping Weight: 31/2 lbs. Master Pack 4 ea.



1400 CASE™

Outside Dimensions: 13" L x 12"W x 6" D. Inside Dimensions: 12"/4" L x 9"/4" W x 5"/4" D. Fold down handle, Integral Personalized Nameplate. Colors: Black, Silver. Shipping Weight: 5"/2 lbs. Master Pack 4 ea.



1450 CASE™

Outside Dimensions: 16" L x 13"W x 7" D. Inside Dimensions: 15" L x 10½"W x 6" D. Fold down handle, Integral Personaliz ile, integral Personalized Nameplate. Colors: Black, Silver. Shipping Weight: 71/2 lbs. Master Pack 4 ea.



1500 CASE™

Outside Dimensions: 181/2" L x 151/4" W x 67/8" D. Inside Dimensions: 17" L x 11 3 /a "W x 6" D. Integral Personalized Nameplate. Colors: Black, Silver. Shipping Weight: 81/2 lbs. Master Pack 4 ea.



1520 CASE™

Outside Dimensions: 191/4" L x 151/2" W x 71/2" D. Inside Dimensions: 18" L x 12 7/8" W x 6 3/4" D.Fold down handle, Large Latch Locks, Integral Personalize Nameplate. Colors: Black, Silver. Shipping Weight: 101/2 lbs. Master Pack 2 ea.

> Technical Safety & Rescue Products

P.O. Box 588 Sandston, VA 23150-0588 Richmond 804-222-1238 800-359-4617




1527 CONVERTIBLE TRAVEL BAG

Designed to fit inside a 1520 case. Converts insi for complete hardknocks travel protection. Just pull out for short day trips etc. sling over shoulder and go without rearranging anything. Made of heavy duty waterproof padded Ballistic Nylon with multiple front and side gadget pockets. Velcro® padded adjustable divider system protects complete outfits. Exclusive top/side access to equipment while in use with thick curved and padded shoulder strap. All nylon zippers and clips. Size: $17^1/2^* \times 12^1/2^* \times x \in 1/4^* D$. Color: Black.




1550 CASE™ 8½" D 20½" Outside Dimensions: 20½" L x 16¾4" W x 8½" D. Inside Dimensions: 19"L x 14"W x 734"D. Fold down handle. Integral Personalized Nameplate. Colors: Black, Silver. Shipping Weight: 14 lbs. Master Pack 2 ea.




1600 CASE™

Outside Dimensions: 231/4" L x 203/4" W x 9" D. Inside Dimensions Nameplate, Colors: Black, Silver. 22" L x 17"W x 8" D. Integral Personalized Shipping Weight: 17 lbs. Master Pack 2 ea.



Outside Dimensions: 31½"L x 20¾" W x 11½" D. Inside Dimensions: 29"L x 18" W x 10½" D. Fold down main handle, large side trunk handles 2" hard rubber transport wheels included, Integral Personalized Nameplate. Colors: Black, Silver. Shipping Weight: 17½ lbs. Master Pack 1 ea.



1700 TRAVEL VAULT II® Outside Dimensions: 37" L x 17½" W x 6" D. Inside Dimensions: 36" L x 13½" W x 5" D. Colors: Black. Shipping Weight: 17 lbs. Master Pack 2 ea

1750 TRAVEL VAULT® 6"D

52" Outside Dimensions: 52" L x 171/2" W x 6" D. Inside Dimensions: 503/4" L x 131/2" W x 5" D. Colors: Black. Shipping Weight: 29 lbs. Master Pack 1 ea.



Unconditional Lifetime Guarantee of Excellence

"Pelican" Cases, Field Support Packaging (Continued)

General Support / Generator

Make / Model:

Description:

The specifications shown are for the generator in use by the WMA MMST. The unit in use is a Dayton heavy duty, long life model powered by a Yanmar air cooled diesel engine. It is rated for continuous use at full load.

Features / Benefits

- Reduces or eliminates reliance on public utilities, allows for independent operation
- Portable, allows for rapid response / deployment

Specifications

- 5kW-rated watts, 6,500-maximum VA
- 120v/41.7 amps, 240v/20.8 amp
- Two 120/240v, 30-amp Twist-Lock outlets
- Two 120V, 20-amp outlets
- 19-HP Yanmar diesel engine
- Electric start
- 4.2-gal. fuel tank / 8.3-hour run time @ 50% duty cycle
- Full wraparound steel frame with handle, lifting eye, and pneumatic wheel kit

Color

N/A

Size

35 in. (L) x 21 in. (W) x 26 in. (H)

Weight

350 lbs.

Shelf Life

N/A

Accessories

Varies with manufacturer / supplier

Approvals

WMA MMST has procured a Dayton HD, 5kW, Diesel powered generator, Model #1N171 from: Grainger Industrial Equipment & Supply Co. (800) 323-0620

Manufacturer / Distributor

Local Purchase

Portable Shelters

Make / Model:

TVI, Rigid-Frame Shelters

Description:

TVI shelters are rapidly deployable, rigid articulating frame shelters with available inner canopy liners and floors. The TVI line of shelters offer side and end entrance doors, and a flexible interior configuration due to the removable compartment dividers. Custom configurations are available and all shelters may be complexed for maximum utility.

Features / Benefits

- All components are replaceable in the field
- No external power sources or tools required for setup / deployment
- Frames are able to withstand up to 30-percent damage and remain self-supporting
- Canopies, liners, and floors are UV resistant and fire retardant

Specifications

- All shelter systems include the shelter, transport and storage bag, and anchor kit
- Various lengths and widths available

Color

Choice of red, green, white, or user may specify

Size

Contact manufacturer for complete list of standard size options

Weight

N/A

Shelf Life

N/A

Accessories

- Specialty shelters
- •Trailer-mounted HVAC systems and generators
- •Mass Casualty Decontamination Systems
- •Contact manufacturer for complete list of accessory / related items

Approvals

N/A

Manufacturer / Distributor

TVI Corporation Glendale, MD 19154 (800) 598-9712

TVI



TVI, Rigid-Frame Shelters, Portable Shelters

Portable Shelter

Make / Model:

Zumro, Air-Beam Inflatable Shelters

Description:

The ZUMRO shelters are portable and rapidly deployable "Airframe" supported shelters with available inner canopy liners and integral floors. The canopies offer multiple entrance doors and a variety of interior configurations. Shelter system components include a shower stall, electric inflator / deflator, two disposable decontamination pools, carrying case, and field repair kit.

Features / Benefits

- Dual inflator connections for quick setup, can be inflated with included inflator, standard air bottles or air compressor
- Canopy and floor replaceable in the field

Specifications

- Air frame and floor material is flame retardant, vinyl coated nylon, canopy material is urethane coated, UV protected nylon, window and skylight material is 30 mil. clear, marine grade vinyl
- Heater / air-conditioner connection sleeve
- Operating pressure is 3 psi, airframe is equipped with pressure relief valve

Color:

White air frame, orange canopy, and black floor

Size

Standard sizes for multi-purpose use; custom sizing available, contact manufacturer for details

Weight

Dependent on size

Shelf Life

N/A

Accessories

N/A


Approvals

N/A

Manufacturer / Distributor

Zumro, Inc. Willow Grove, PA 19090 (800) 932-6003 RMC Medical Philadelphia, PA 19154 (215) 824-4100

Zumro Portable Tent / Shelter



H66-13

Zumro, Air-Beam Inflatable Shelters

Portable Shelters

Make / Model:

Air-Lite, Air-Beam Inflatable Shelters

Description:

The Air-Lite shelters are rapidly deployable "Airframe" supported shelters with double doors, integrated screen doors, and integrated floors. The shelters may be complexed to create custom configurations.

Features / Benefits

- Heavyweight urethane coated nylon fabric construction
- Multiple inflator options for quick setup, can be inflated with included inflator, standard air bottles, or air compressor
- Floor replaceable in the field

Specifications

- Three sizes available:
 - 12 ft. x 15 ft. w/ 5 ft. sidewall, 88 lbs. (approx.)
 - 10 ft. x 15 ft. w/ 4 ft. sidewall, 75 lbs. (approx.)
 - 10 ft. x 10 ft. w/ 4 ft. sidewall, 60 lbs. (approx.)

Color

Blue and Silver

Size

N/A

Weight

See above

Shelf Life

N/A

Accessories

- Wheeled storage box
- 12-volt de Electric Inflator
- 12-volt Portable Battery
- 110-volt ac electric Inflator

Approvals

N/A

Manufacturer / Distributor

Air-Lite Inflatables, Inc. Haines, OR 19154 (800) 200-1960





Air-Lite, Air-Beam Inflatable Shelters